



**nsMCDLibrary**  
**Neuroshare implementation for MC\_Rack data**

---

Multi Channel Systems MCS GmbH  
Aspenhastrasse 21  
72770 Reutlingen  
Germany  
Fon +49-71 21-90 92 5 – 0  
Fax +49-71 21-90 92 5 – 11  
[info@multichannelsystems.com](mailto:info@multichannelsystems.com)  
[www.multichannelsystems.com](http://www.multichannelsystems.com)

nsMCDLibrary Version: 2.7  
Stand: 2009-11-17  
Autor: Jens Pätzold

## ***General***

To access MC\_Rack data files from the Neuroshare API, the data are mapped to the structure given by Neuroshare API. This API is implemented in a Windows DLL or Linux shared library by a set of predefined functions.

Up to now, not all possible data streams within a MC\_Rack file could be read by the Neuroshare API.

The nsMCDLibrary.dll is implemented according to the Neuroshare API Version 1.3.

This paper only handles what is specific to data from MC\_Rack. For the documentation of the Neuroshare API see there: <http://neuroshare.sourceforge.net/>

## ***Installation***

The nsMCDLibrary.dll can be downloaded from our webpage at <http://www.multichannelsystems.com/downloads>. For Windows the nsMCDLibrary\_2.7.zip file contains the nsMCDLibrary.dll itself, this documentation and examples for FIND (<http://find.bccn.uni-freiburg.de/>).

Just unzip the file anywhere on your computer. If your FIND installation has only an older version of the nsMCDLibrary.dll, just replace the nsMCDLibrary.dll with the newer one.

For Linux both 32bit and 64bit versions are available in the tar file nsMCDLibrary\_Linux32\_2.7.tar.gz or nsMCDLibrary\_Linux64\_2.7.tar.gz. They contain the shared library, a small c example file and this documentation.

For MacOSX a library (nsMCDLibrary\_MacOSX.tar.gz) exists with the small c example file and this documentation.

## ***Implementation details***

MC\_Rack has two different recording types. Both continuous and triggered data can be handled with the neuroshare library. The mapping is a little bit different for each type, so the mapping is shown for each type separately.

### **Continuous data mapping**

<b>Stream</b>	<b>Short name</b>	<b>Entity</b>
Electrode Raw Data	elec	Analog Entity
Analog Raw Data	anlg	Analog Entity
Filtered Data	filt	Analog Entity
Channel Tool Data	chtl	Analog Entity
Digital Data	digi	Analog Entity, each binary digit as Event Entity
Spikes	spks	Segment Entity
Trigger	trig	Event Entity

Other stream are not mapped up to now.

### **Triggered data mapping**

<b>Stream</b>	<b>Short name</b>	<b>Entity</b>
Electrode Raw Data	elec	Analog Entity and Segment Entity
Analog Raw Data	anlg	Analog Entity and Segment Entity
Filtered Data	filt	Analog Entity and Segment Entity
Channel Tool Data	chtl	Analog Entity and Segment Entity
Digital Data	digi	Analog Entity and Segment Entity, each binary digit as Event Entity
Spikes	spks	Segment Entity
Trigger	trig	Event Entity

Please be aware that some streams are mapped to two entities. The contents is essentially the same. Other stream are not mapped up to now.

### **Missing streams**

All other possible streams as average and different parameters are up to now not implemented in the nsMCSLibrary.dll.

## ***Entities***

### **Event entities**

Each of the 16 bits of the digital data channel is represented as one event entity. Each change on the individual digital line results on a new event in the entity with the next index. The event information is in one 8bit (byte) value. Its value is 1 for a change from 0 to 1 on the digital line and -1 (255) for a change from 1 to 0 on the digital line.

Triggers are mapped to one event entity for each trigger stream. The event is represented in two 16bit (word) values. The first one gives the slope of the Trigger, the second the Trigger value.

### **Analog entities**

Raw electrode data, analog data and filtered data are mapped to analog entities. Each channel gives one analog entity in Neuroshare. The different gain is considered. The data are given in Volt of the original input.

The 16 bit of the digital data channel are represented in one analog entity as raw data with a range from 0 to 65535, with minimum step size of one. Each sample value is included in the entity.

For continuous data the first index is starting at time 0. Then all data are following with the next index for the next sample point. The time difference of each index is the reciprocal value of the sample rate.

For triggered data the indices are not always continuous in time. If you are reading many indices at one time, you get, according to the Neuroshare specifications, in `pdwContCount` the number of samples that are continuous:

„Although the samples in an analog entity are indexed, they are not guaranteed to be continuous in time and may contain gaps between some of the indexes. When the requested data is returned, `pdwContCount` contains the number of Analog items, starting from `dwStartIndex`, which do not contain a time gap“. [NeuroshareAPI-1-3.pdf]

For analog entities there is no timestamp directly included in the API for reading the data. You get the timestamp with the function `ns_GetTimeByIndex`.

### **Segment entities**

Spike streams are represented in segment entities. Spikes (or wave forms) are sampled independently for each channel. This means there is always only one source per segment entity. The number of sources is always 1. Each original channel from MC\_Rack data forms one entity in Neuroshare. The different gain is considered. The data are given in Volt of the original input.

For triggered data all analog entities are also represented as segment entities. One segment contains the data of one sweep for one channel.

### **Neural event entities**

There are no streams in MC\_Rack data that are mapped to neural events.

## ***Naming convention for entities***

The entities are named as follow:

- 8 letter: stream name with:
  - 4 letters name
  - 4 digits stream number.
- Space
- 4 letters: HWID
- Space
- 4 digits: index of channel
- Space
- 8 letters: name of channel (right justified)

From here on, this applies only to event entities of digital data:

- Space
- 2 digits sub channel for events of digital data

for example:

```
aaaa0000 0000 0000 xxxxxxxx 00
elec0001 0001 0000          21
digi0001 0063 0001          D1
digi0001 0063 0001          D1 02
```

## ***Missing***

The function `ns_GetLastErrorMsg` is not implemented yet.