

## # \$ K + The ABF Header Section

The ABF header is the first block of data at the start of an ABF data [file](#). The header contains parameters that describe the stimulation, the acquisition and the [hierarchy](#) of the data. It describes the contents of the data file and contains entries to describe the settings in effect when the data file was acquired.

In version 1.8, the header is 6144 bytes long. This will most likely be increased in subsequent versions. Third party programs should NOT rely on the size of the header, or retrieve information based on a byte offset. Only use the documented variables defined in the header.

See the file [ABFHEADR.H](#) for a "C" definition of the ABFFileHeader structure.

### ABF Header Sub Sections

- |  |   |
|--|---|
| 1. <a href="#">File ID and Size information</a>  | 17. <a href="#">Extended File Structure</a>                       |
| 2. <a href="#">File Structure</a>                | 18. <a href="#">Extended Multi-Channel Information</a>            |
| 3. <a href="#">Trial Hierarchy Information</a>   | 19. <a href="#">Train Parameters</a>                              |
| 4. <a href="#">Display Parameters</a>            | 20. <a href="#">Extended Epoch Waveform and Pulses</a>            |
| 5. <a href="#">Hardware Information</a>          | 21. <a href="#">Extended DAC Output File</a>                      |
| 6. <a href="#">Environmental information</a>     | 22. <a href="#">Extended Pre-sweep (Conditioning) Pulse Train</a> |
| 7. <a href="#">Multi-channel information</a>     | 23. <a href="#">Extended Variable Parameter User List</a>         |
| 8. <a href="#">Synchronous timer outputs</a>     | 24. <a href="#">Extended On-Line Subtraction</a>                  |
| 9. <a href="#">Epoch Waveform and Pulses</a>     | 25. <a href="#">Extended Environmental Information</a>            |
| 10. <a href="#">DAC Output File</a>              | 26. <a href="#">Extended Statistics Measurements</a>              |
| 11. <a href="#">Conditioning pulse train</a>     | 27. <a href="#">Application Version Data</a>                      |
| 12. <a href="#">Variable Parameter User List</a> | 28. <a href="#">LTP Protocol</a>                                  |
| 13. <a href="#">Autopeak measurement</a>         | 29. <a href="#">Digidata 132x Trigger Out</a>                     |
| 14. <a href="#">Channel Arithmetic</a>           | 30. <a href="#">Epoch Resistance</a>                              |
| 15. <a href="#">On-line Subtraction</a>          | 31. <a href="#">Alternating Episodic Mode</a>                     |
| 16. <a href="#">Miscellaneous Parameters</a>     | 32. <a href="#">Post-Processing Actions</a>                       |

The total header length = 6144 bytes.

Offsets in the lists referenced above are Byte offsets.

### Notes:

[Data Types](#)

[Version Numbers](#)

[ADC Channel Numbering](#)

[Composite Values](#)

[Indexing Arrays in the ABF Header](#)

---

# The\_ABF\_Header\_Section

\$ The ABF Header Section

K The ABF Header Section

+ ABFFUNC

## # \$ K + ADC Channel Numbering

Axon's data acquisition programs distinguish between physical and logical channel numbers. Physical channel numbers are the channel numbers used internally to communicate with the acquisition hardware. Logical channel numbers are the external connector labels on the front panel of the acquisition hardware. Logical channel numbers are used only for presentation to the user. Physical channel numbers are used everywhere else. For example, parameters are stored using physical channel number order (0 to 15) for such structures as the sampling sequence array and the entries for the external lowpass and highpass filters. Similarly, a physical channel number is used for the Trigger channel. At the time of printing this document, the only digitizer known to have different physical and logical channel numbering is the (now obsolete) TL-2 interface.

---

# ADC\_Channel\_Numbering

\$ ADC Channel Numbering

K ADC Channel Numbering

+ ABFFUNC

## # \$ K + Indexing Arrays in the ABF Header

To get a Logical channel number from a Physical channel number, simply index the nADCPtoLchannelMap array by the Channel number you wish to convert. Thus nADCPtoLchannelMap[1] provides the Logical Channel Number for Physical Channel Number 1. This array is always symmetrical, so it can be used in the same way to convert back to Physical Channel Numbers from Logical Channel Numbers.

The first thing to look at is the nADCSamplingSeq array. This tells you which physical [ADC](#) channels were acquired and in what order. The first entry in this array is the Physical channel number of the first ADC channel acquired, followed by the second etc. There are nADCNumChannels channels in this array. All ADC arrays except for the nADCSamplingSeq are indexed through Physical channel numbers. These include: sADCChannelName, sADCUnits, etc.

Special Note: All array indexing within the header and within the ABF routines start at 0, except Sweep number, which starts at 1.

```
void ShowFirstAcquiredChannelInfo( int nFile, ABFFileHeader *pFH )
{
    // Get first physical channel number and name
    int nFirstPhysicalChannel = pFH->nADCSamplingSeq[0];
    char *psSignalName = pFH->sADCChannelName[nFirstPhysicalChannel];

    // Get channel number to show to user
    int nFirstLogicalChannel = pFH->nADCPtoLchannelMap[nFirstPhysicalChannel];
    fprintf("The first acquired channel (%s) comes from ADC channel %d\n",
        psSignalName, nFirstLogicalChannel);

    // Get start time of first episode of first channel
    int nFirstEpisode = 1;
    float fStartTime;
    ABF_GetStartTime(nFile, pFH, nFirstPhysicalChannel, nFirstEpisode,
        &fStartTime, NULL);
}
```

---

# Indexing\_Arrays\_in\_the\_ABF\_Header

\$ Indexing Arrays in the ABF Header

K Indexing Arrays in the ABF Header

+ ABFHELP

## # \$ K + Composite Values

Axon Instruments' acquisition programs can automatically read the gain and lowpass filter settings of the analog [instrument](#) (e.g. patch clamp) on upto 16 [ADC](#) channels. If the ABF reading routines are used the telegraphed gain must be taken into account when calculating the signal scale factor for the telegraphed channel, as must the telegraphed lowpass filter setting. Both of these factors are automatically taken into account by the ABF reading routines.

In some cases, several parameters must be considered together to correctly interpret the data.

Composite ADC scale factor  $f_{\text{InstrumentScaleFactor}}[n] * f_{\text{TelegraphAdditGain}} * f_{\text{ADCProgrammableGain}}[n] * f_{\text{SignalGain}}[n]$ .

Note, the fTelegraphAdditGain only applies to the telegraphed ADC channel.

Composite ADC offset       $f_{\text{InstrumentOffset}}[n] + f_{\text{SignalOffset}}[n]$

Composite lowpass filter -3 dB frequency	A complex function of fSignalLowpassFilter and fTelegraphFilter, for the telegraphed ADC channel. On other ADC channels, only the fSignalLowpassFilter contributes to the composite value.
--	--

The functions [ABFH\\_GetADCtoUUFactors\(\)](#) and [ABFH\\_GetDACtoUUFactors\(\)](#) should be used to get the combined scale factors required for converting ADC/DAC values to [UserUnits](#) and vice versa.

### # Composite\_Values

<sup>\$</sup> Composite Values

<sup>K</sup> Composite Values<sup>+</sup> ABFFUNC

## # \$ K + **Data Types**

Unused real parameters should be filled with zeros. Unused strings should be filled with the space character (ASCII #32).

Parameters for unsampled ADC channels should be filled with the indicated default.

---

# Data\_Types

\$ Data Types

K Data Types

+ ABFFUNC

## # \$ K + **Version Numbers**

The file version number consists of a major and a minor number. For example, the "1" in Version 1.0 is the major number, and the "0" is the minor number.

In general, the major version number is updated when changes affect the byte offset of the existing parameters or would otherwise make the file unusable by existing programs. The minor version number is updated when unused parameters are utilized. In most cases, existing programs will not be affected since they should not be dependent upon the unused parameters.

---

# Version\_Numbers

\$ Version Numbers

K Version Numbers

+ ABFFUNC

## # \$ K + **ABF Hardware and Storage Limits**

Some of the more important limitations in the range of hardware supported and the size of components in ABF formatted [files](#) are listed here.

- Sixteen physical [ADC](#) channels, numbered 0-15.
- Up to sixteen bits per ADC word.
- Up to 1,032,258 multiplexed samples per [sweep](#) in high-speed oscilloscope mode, [fixed-length](#) event-driven mode and [episodic stimulation](#) mode.
- Up to 2 G multiplexed samples per segment in [variable-length](#) event-driven mode and [gap-free](#) mode.
- Stimulus waveform can be generated on two [DACs](#) channel simultaneously.
- Pre-sweep train (was previously called conditioning train) can be generated on two [DACs](#) channel simultaneously.
- One Math Channel.
- Upto 16 telegraphed [instruments](#)
- Leak subtraction can be applied to ADC channels simultaneously.
- One set of display amplifications and offsets.
- Only one averaged Run per file!

---

# ABF\_Hardware\_and\_Storage\_Limits

\$ ABF Hardware and Storage Limits

K ABF Hardware and Storage Limits

+ ABFHELP

**Special Note:** A header entry name prefixed with underscore "\_" indicates that this parameter is discontinued or has been replaced with an updated version. In such cases the parameter is also tagged as **Discontinued** in the description field and where appropriate a link has been added to the updated version of the parameter in the Word version of this document.

## # \$ K + File ID and Size information

(40 bytes).

These ten entries represent two types of information: (1) File identification information; (2) Parameters whose actual value may be different to the value requested before the acquisition commenced, e.g. due to a user abort. The requested values are located elsewhere.

Offset	Header Entry Name	Type	Description
0	IFileSignature	<a href="#">long</a>	File type used for format identification. Possible values are: "ABF ", "CLPX" and "FTCX". Used to create the numbers in nFileType. (In old pCLAMP and Axotape data files, the first four bytes were a <a href="#">float</a> : 1 = CLAMPEX, 10 = FETCHEX/AxoTape. This is translated on reading into either CLPX or FTCX as appropriate.)
4	fFileVersionNumber	<a href="#">float</a>	File format version stored in the data file during acquisition. The present version is 1.5. (In old pCLAMP and Axotape data files, this parameter is in the range 2.0-5.3)
8	nOperationMode	<a href="#">short</a>	Operation mode: 1 = Event-driven, variable length; 2 = Oscilloscope, loss free (Same as Event-driven, <a href="#">fixed length</a> ); 3 = Gap-free; 4 = Oscilloscope, high-speed; 5 = episodic stimulation (Clampex only).
10	IActualAcqLength	<a href="#">long</a>	Actual number of <a href="#">ADC</a> samples (aggregate) in data file. See IAcqLength. Averaged <a href="#">sweeps</a> are included.
14	nNumPointsIgnored	<a href="#">short</a>	Number of points ignored at data start. Normally zero, but non-zero for <a href="#">gap-free</a> acquisition using AXOLAB configurations with one or more ADS boards.
16	IActualEpisodes	<a href="#">long</a>	Actual number of <a href="#">sweeps</a> in the file including averaged sweeps. See IEpisodesPerRun. If nOperationMode = 3 ( <a href="#">gap-free</a> ) the value of this parameter is 1.
20	IFileStartDate	<a href="#">long</a>	Date when data portion of this file was first written to. Stored as YYMMDD. If YY is in the range 80-99, prefix with "19" to get the year. If YY is in the range 00-79, prefix with "20" to get the year.
24	IFileStartTime	<a href="#">long</a>	Time of day in seconds past midnight when data portion of this file was first written to.

---

# File\_ID\_and\_Size\_information

\$ File ID and Size information

K File ID and Size

information;IFileSignature;fFileVersionNumber;nOperationMode;IActualAcqLength;nNumPointsIgnored;IActualEpisodes;IFileStartDate;IFileStartTime;IStopwatchTime;fHeaderVersionNumber;nFileType;nMSBinFormat

+ ABFFUNC



28	lStopwatchTime	<u>long</u>	Time since the stopwatch was zeroed that the data portion of this file was first written to. Not supported by all programs. Default = 0.
32	fHeaderVersionNumber	<u>float</u>	Version number of the header structure returned by the ABF_ReadOpen function. Currently 1.8. This parameter does not identify the data file format. See fFileVersionNumber.
36	nFileType	<u>short</u>	Numeric equivalent of file type. 1 = ABF file; 2 = Old FETCHEX file (FTCX); 3 = Old Clampex file (CLPX). See sFileType.
38	nMSBinFormat	<u>short</u>	Storage method for real numbers in the header. Also see nDataFormat. 0 = IEEE format; 1 = Microsoft Binary format (old files only).

## # \$ K + File Structure

(78 bytes).

Header entries describing the structure of the [file](#). See Extended File Structure section for new variables

Offset	Header Entry Name	Type	Description
40	IDataSectionPtr	<a href="#">long</a>	Block number of start of Data section.
44	ITagSectionPtr	<a href="#">long</a>	Block number of start of Tag section.
48	INumTagEntries	<a href="#">long</a>	Number of Tag entries.
52	IScopeConfigPtr	<a href="#">long</a>	Block number of the ABF Scope Config section. (was block number of start of Long Description section.)
56	INumScopes	<a href="#">long</a>	Number of ABFScopeConfig structures in the ABF Scope Config section. (was number of lines of Long Description.)
60	_IDACFilePtr	<a href="#">long</a>	Block number of start of <a href="#">DAC</a> file section. No longer used.  <b>Discontinued:</b> Use <a href="#">IDACFilePtr</a>
64	_IDACFileNumEpisodes	<a href="#">long</a>	Number of <a href="#">sweeps</a> in the <a href="#">DAC</a> file section. Sweeps are not multiplexed.  <b>Discontinued:</b> Use <a href="#">IDACFileNumEpisodes</a>
68	SUnused001	4char	Unused.
72	IDeltaArrayPtr	<a href="#">long</a>	Block number of start of Delta Array section.
76	INumDeltas	<a href="#">long</a>	Number of entries in Delta Array section.
80	IVoiceTagPtr	<a href="#">long</a>	Block number of start of Voice Tag section.
84	IVoiceTagEntries	<a href="#">long</a>	Number of Voice Tag entries.
88	IUnused002	<a href="#">long</a>	(was number of automatic entries in Notebook section)
92	ISynchArrayPtr	<a href="#">long</a>	Block number of start of the Synch Array section.
96	ISynchArraySize	<a href="#">long</a>	Number of pairs of entries in the Synch Array section. If averaging is enabled, this includes the entry for the averaged <a href="#">sweep</a> .
100	nDataFormat	<a href="#">short</a>	Data representation. 0 = 2-byte integer; 1 = IEEE 4 byte <a href="#">float</a> .
102	nSimultaneousScan	<a href="#">short</a>	ADC Channel Scanning Mode: 0=Multiplexed; 1=Simultaneous Scanning (currently unimplemented)
104	IStatisticsConfigPtr	<a href="#">long</a>	Block number of start of Scope Config section.
108	IAnnotationSectionPtr	<a href="#">long</a>	Block number of start of annotations section
112	INumAnnotations	<a href="#">long</a>	Number of annotations
116	sUnused003	<a href="#">2char</a>	Unused.

---

# File\_Structure

\$ File Structure

K File

Structure;IDataSectionPtr;ITagSectionPtr;INumTagEntries;IScopeConfigPtr;INumScopes;IDACFilePtr;IDACFileNumEpisodes;IDeltaArrayPtr;INumDeltas;ISynchArrayPtr;ISynchArraySize;nDataFormat;

+ ABFFUNC

## # \$ K + Trial Hierarchy Information

(82 bytes).

Header entries describing the [trial hierarchy](#).

Offset	Header Entry Name	Type	Description
118	channel_count_acquired	<a href="#">short</a>	Number of analog input channels acquired. Can be less than nADCNumChannels. Currently unimplemented.
120	nADCNumChannels	<a href="#">short</a>	Number of analog input channels sampled.
122	fADCSampleInterval	<a href="#">float</a>	Interval between multiplexed <a href="#">A/D</a> samples (us). In Clampex, this is also known as the first clock interval.
126	fADCSecondSampleInterval	<a href="#">float</a>	Second definition of the interval between multiplexed A/D samples (us). If this interval is zero, the first clock interval is used for the whole of each <a href="#">sweep</a> . If this interval is non-zero, the second clock interval starts at the sample number specified by IClockChange.
130	fSynchTimeUnit	<a href="#">float</a>	Time unit for start time in the Synch Array section: 0 = Value in sample intervals; nn = Value in $\mu$ s. See notes below for the Synch Array section. pCLAMP 6.0 and AxoTape 2.0 use value in sample intervals only.
134	fSecondsPerRun	<a href="#">float</a>	Requested acquisition length in seconds. 0 means use available disk space. -1 means ignore this parameter and refer to IEpisodesPerRun.
138	INumSamplesPerEpisode	<a href="#">long</a>	Number of multiplexed <a href="#">ADC</a> samples per <a href="#">sweep</a> if nOperationMode is 2, 4, or 5. Undefined if nOperationMode is 1 or 3.
142	IPreTriggerSamples	<a href="#">long</a>	Pre-trigger interval stored as number of multiplexed ADC samples (note that this is the underlying number of ADC samples, not the number of pre-trigger samples in the trigger channel). FETCHEX uses the same value for the post-trigger interval. Undefined if nOperationMode is 3 or 5.
146	IEpisodesPerRun	<a href="#">long</a>	Requested number of sweeps/ <a href="#">run</a> . 0=Run until terminated by user. If nOperationMode = 3 ( <a href="#">gap free</a> ), this parameter is 1 and the requested acquisition length is set in fSecondsPerRun.
150	IRunsPerTrial	<a href="#">long</a>	Requested number of runs/ <a href="#">trial</a> . 0=Run until terminated by user. Runs are averaged. If nOperationMode = 3 (gap free), the value of this parameter is 1. See IAverageCount.
154	INumberOfTrials	<a href="#">long</a>	Number of trials to be acquired. Note, only one trial is contained in a data file. This parameter is used for Clampex acquisition control. -1 = continuous.

---

# Trial\_Hierarchy\_Information

\$ Trial Hierarchy Information

K Trial Hierarchy

Information;nADCNumChannels;fADCSampleInterval;fADCSecondSampleInterval;fSynchTimeUnit;fSecondsPerRun;INumSamplesPerEpisode;IPreTriggerSamples;IEpisodesPerRun;IRunsPerTrial;INumberOfTrials;nAveragingMode;nUndoRunCount;nFirstEpisodeInRun;fTriggerThreshold;nTriggerSource;nTriggerAction;nTriggerPolarity;fScopeOutputInterval;fEpisodeStartToStart;fRunStartToStart;fTrialStartToStart;IAverageCount;IClockChange;nAutoTriggerStrategy

+ ABFFUNC

158	nAveragingMode	<a href="#">short</a>	Data averaging strategy: 0 = No averaging; 1 = Average but don't save raw data; 2 = Average and save raw data.
160	nUndoRunCount	<a href="#">short</a>	Frequency with which temporary undo file created during averaging is updated. -1 = Disabled; 0 = after each pseudo-doubling; N = after every N runs.
162	nFirstEpisodeInRun	<a href="#">short</a>	First <a href="#">sweep</a> number (normally 1).
164	fTriggerThreshold	<a href="#">float</a>	Trigger threshold value for Event or Oscilloscope mode (user units).
168	nTriggerSource	<a href="#">short</a>	Trigger source: N ( $\geq 0$ ) = Physical channel number selected for threshold detection; -1 = external trigger; -2 = keyboard; -3 = use start-to-start interval.  If nOperationMode=3 (gap-free) 0= start immediately.
170	nTriggerAction	<a href="#">short</a>	Trigger start strategy: 0 = start one <a href="#">sweep</a> ; 1 = start one <a href="#">run</a> ; 2 = start one <a href="#">trial</a> .
172	nTriggerPolarity	<a href="#">short</a>	Trigger on: 0 = rising edge of waveform; 1 = falling edge. The terms "rising" and "falling" refer to the direction of the data transitions on the screen during acquisition.
174	fScopeOutputInterval	<a href="#">float</a>	Interval between digital output used for external oscilloscope, or equivalent, and the start of each <a href="#">sweep</a> (ms)
178	fEpisodeStartToStart	<a href="#">float</a>	Time between start of sweeps (seconds). Use when nTriggerSource = "start-to-start".
182	fRunStartToStart	<a href="#">float</a>	Time between start of runs (seconds).
186	fTrialStartToStart	<a href="#">float</a>	Time between trials (seconds).
190	IAverageCount	<a href="#">long</a>	The actual number of runs contributing to each average. See IEpisodesPerRun or IRunsPerTrial.
194	IClockChange	<a href="#">long</a>	The multiplexed ADC sample number after which the second sampling interval commences. 0 = INumSamplesPerEpisode/2. See fADCSecondSampleInterval.
198	nAutoTriggerStrategy	<a href="#">short</a>	0 = Do not auto trigger; 1 = Autotrigger. (Only significant when nOperationMode == ABF_HIGHSPEEDOSC)

## # \$ K + Display Parameters

(44 bytes).

Header entries describing display options in effect during acquisition.

Offset	Header Entry Name	Type	Description
200	nDrawingStrategy	<a href="#">short</a>	Strategy for the drawing of raw data: 0 = not at all; 1 = update immediately as data is acquired; 2 = update at the end of each <a href="#">sweep</a> or <a href="#">trace</a> ; 3 = update at the end of each <a href="#">run</a> ;
202	nTiledDisplay	<a href="#">short</a>	<a href="#">ADC</a> channel display arrangement: 0 = Superimposed; 1 = Tiled.
204	nEraseStrategy	<a href="#">short</a>	Automatically erase screen: 0 = not at all; 1 = before each sweep; 2 = before each run.  <b>Discontinued:</b> Use scope config entry instead.
206	nDataDisplayMode	<a href="#">short</a>	Data display mode: 0 = Points; 1 = Lines.
208	IDisplayAverageUpdate	<a href="#">long</a>	Display averaged data: -1 = at end of <a href="#">trial</a> ; 0 = after each pseudo-doubling; N = after N runs.
212	nChannelStatsStrategy	<a href="#">short</a>	Show channel statistics in <a href="#">gap-free</a> mode: 0 = No; 1 = Yes.
214	ICalculationPeriod	<a href="#">long</a>	Length of the real-time statistics update period in samples. Conventionally a multiple of 1024. Default = 16384.  <b>Discontinued:</b> Use <a href="#">fStatisticsPeriod</a>
218	ISamplesPerTrace	<a href="#">long</a>	Number of multiplexed ADC samples in displayed <a href="#">trace</a> .
222	IStartDisplayNum	<a href="#">long</a>	Starting sample number for <a href="#">sweep</a> display: N = starting sample number. (Use N = 1 to start from the first sample.)
226	IFinishDisplayNum	<a href="#">long</a>	Finishing sample number for sweep display: 0 = finish at end of sweep; N = finishing sample number.
230	nMultiColor	<a href="#">short</a>	Color control for multi-trace displays. 0 = single color for all traces; 1 = two or more colors for traces.
232	nShowPNRawData	<a href="#">short</a>	0 = display corrected data; 1 = display raw data.
234	fStatisticsPeriod	<a href="#">float</a>	Length of the real-time statistics update period in seconds.  Not used.
238	IStatisticsMeasurements	<a href="#">long</a>	Bit mask for statistics measurements to display. This applies to trigger statistics only. For episodic statistics use IStatsMeasurements

Above Threshold: 0x00000001

Event Frequency: 0x00000002

---

# Display\_Parameters

\$ Display Parameters

K Display

Parameters;nDrawingStrategy;nTiledDisplay;nEraseStrategy;nDataDisplayMode;IDisplayAverageUpdate;nChannelStatsStrategy;ICalculationPeriod;ISamplesPerTrace;IStartDisplayNum;IFinishDisplayNum;nMultiColor;nShowPNRawData;fStatisticsPeriod

+ ABFFUNC

Mean Open Time: 0x00000004

Mean Closed Time: 0x00000008

242

nStatisticsSaveStrategy

short

Strategy used to save statistics: No Auto Save = 0; Auto  
Save = 1

## # \$ K + Hardware Information

(16 bytes).

Description of hardware properties of digitizer used for acquisition.

Offset	Header Entry Name	Type	Description
244	fADCRange	<a href="#">float</a>	<a href="#">ADC</a> positive full-scale input in volts (e.g. 10.00V).
248	fDACRange	<a href="#">float</a>	<a href="#">DAC</a> positive full-scale range in volts.
252	IADCResolution	<a href="#">long</a>	Number of <a href="#">ADC</a> counts corresponding to the positive full-scale voltage in ADCRange (e.g. 2000, 2048, 32000 or 32768).
256	IDACResolution	<a href="#">long</a>	Number of <a href="#">DAC</a> counts corresponding to the positive full-scale voltage in DACRange.

---

# Hardware\_Information

\$ Hardware Information

K Hardware Information;fADCRange;fDACRange;IADCResolution;IDACResolution

+ ABFFUNC

## # \$ K + Environmental information

(118 bytes).

Description of telegraph properties and comment info. See Extended Environmental Information for new parameters.

Offset	Header Entry Name	Type	Description
260	nExperimentType	<a href="#">short</a>	Experiment type: 0 = Voltage Clamp; 1 = Current Clamp.
262	_nAutosampleEnable	<a href="#">short</a>	Enable storage of autosample information: 0 = Disabled; 1 = Automatic; 2 = Manual.  <b>Discontinued:</b> Use <a href="#">nTelegraphEnable</a>
264	_nAutosampleADCNum	<a href="#">short</a>	Physical <a href="#">ADC</a> channel number to which autosampled parameters apply.  <b>Discontinued.</b>
266	_nAutosampleInstrument	<a href="#">short</a>	Autosample <a href="#">instrument</a> : (IHS-1 telegraphs are not supported.) Note, for most programs this is an information-only field. For example, in Clampex the autosample instrument is chosen as a configuration item and copied into this field.  <b>Discontinued:</b> Use <a href="#">nTelegraphInstrument</a>
268	_fAutosampleAdditGain	<a href="#">float</a>	Additional gain multiplier of Instrument connected to nAutosampleADCNum. (Optionally autosampled by some acquisition programs.) (Default = 1)  <b>Discontinued:</b> Use <a href="#">fTelegraphAdditGain</a>
272	_fAutosampleFilter	<a href="#">float</a>	Lowpass filter cutoff frequency of Instrument connected to nAutosampleADCNum. (Optionally autosampled by some acquisition programs.) (Default = 100000)  <b>Discontinued:</b> Use <a href="#">fTelegraphFilter</a>
276	_fAutosampleMembraneCap	<a href="#">float</a>	Patch-clamp membrane capacitance compensation. (Optionally autosampled by some acquisition programs.)  <b>Discontinued:</b> Use <a href="#">fTelegraphMembraneCap</a>
280	nManualInfoStrategy	<a href="#">short</a>	Strategy for writing the manually entered information: 0 = Do not write; 1 = Write each <a href="#">trial</a> ; 2 = Prompt each trial..
282	fCellID1	<a href="#">float</a>	Numeric identifier #1, e.g. cell identifier.
286	fCellID2	<a href="#">float</a>	Numeric identifier #2, e.g. temperature in °C.
290	fCellID3	<a href="#">float</a>	Numeric identifier #3.
294	sCreatorInfo	16 <a href="#">char</a>	Name and version of program used to create the file. For example, "AxoTape 2.0" or "Clampex 6.0".

---

# Environmental\_information

\$ Environmental information

K Environmental

information;nExperimentType;nAutosampleEnable;nAutosampleADCNum;nAutosampleInstrument;fAutosampleAdditGain;fAutosampleFilter;fAutosampleMembraneCap;nManualInfoStrategy;fCellID1;fCellID2;fCellID3;sCreatorInfo;sFileComment

+ ABFFUNC



310	_sFileComment	56 <u>char</u>	56 byte ASCII comment string. <b>Discontinued:</b> Use <a href="#">sFileComment</a>
366	nFileStartMillisecs	<u>short</u>	Milliseconds portion of IFileStartTime
368	nCommentsEnable	<u>short</u>	Enable comments field
370	sUnused003a	8 <u>char</u>	Unused.

## # \$ K + Multi-channel information

(1044 bytes).

Description of ADC channel properties.

Offset	Header Entry Name	Type	Description
378	nADCPtoLChannelMap(0-15)	<a href="#">short</a>	<a href="#">ADC</a> physical-to-logical channel map. The entries are in the physical order 0, 1, 2, ... 14, 15. If there are fewer than 16 logical channels in the system, the array is padded with -1. All channels supported by the hardware are present, even if only a subset is used. For example, for the TL-2 the entries would be 7, 6, 5, 4, 3, 2, 1, 0, -1 ... -1.
410	nADCSamplingSeq(0-15)	<a href="#">short</a>	ADC channel sampling <a href="#">sequence</a> . This is the order in which the physical ADC channels are sampled. If fewer than the maximum number of channels are sampled, pad with -1. For example, if two channels are sampled on the TL-2, this array will contain 6, 7, -1 ... -1. If two channels are sampled on the TL-1, this array will contain 14, 15, -1 ... -1.
442	sADCChannelName(0-15)	10char	ADC channel name in physical channel number order. Default = spaces.
602	sADCUnits(0-15)	8char	The user units for ADC channels in physical channel number order. Default = spaces.
730	fADCProgrammableGain(0-15)	<a href="#">float</a>	ADC programmable gain in physical channel number order (dimensionless). Default = 1.
794	fADCDisplayAmplification(0-15)	<a href="#">float</a>	ADC channel display amplification in physical channel number order (dimensionless.) Default = 1.
858	fADCDisplayOffset(0-15)	<a href="#">float</a>	ADC channel display offset in physical channel number order (user units). Default = 0.
922	fInstrumentScaleFactor(0-15)	<a href="#">float</a>	Instrument scale factor in physical ADC channel number order (Volts at ADC / user unit). (Programs would normally display this information to the user as user units / volt at ADC).
986	fInstrumentOffset(0-15)	<a href="#">float</a>	Instrument offset in physical ADC channel number order (user units corresponding to 0 V at the ADC). Default is zero.
1050	fSignalGain(0-15)	<a href="#">float</a>	Signal conditioner gain in physical ADC channel number order (dimensionless). Default = 1.
1114	fSignalOffset(0-15)	<a href="#">float</a>	Signal conditioner offset in physical ADC channel number order (user units). Default = 0.
1178	fSignalLowpassFilter(0-15)	<a href="#">float</a>	Signal-conditioner lowpass filter corner frequency in physical ADC channel number order (Hz). 100000 means lowpass filter is bypassed (i.e. wideband). Default = 100000.

---

# Multi\_channel\_information

\$ Multi-channel information

K Multi-channel

information;nADCPtoLChannelMap;nADCSamplingSeq;sADCChannelName;sADCUnits;fADCProgrammableGain;fADCDisplayAmplification;fADCDisplayOffset;fInstrumentScaleFactor;fInstrumentOffset;fSignalGain;fSignalOffset;fSignalLowpassFilter;fSignalHighpassFilter;sDACChannelName;sDACChannelUnits;fDACScaleFactor;fDACHoldingLevel;nSignalType

+ ABFFUNC

1242	fSignalHighpassFilter(0-15)	<u>float</u>	Signal-conditioner highpass filter corner frequency in physical ADC channel number order (Hz). 0 means highpass filter is bypassed (i.e. DC coupled). -1 means inputs are grounded. Default = 0.
1306	sDACChannelName(0-3)	10 <u>char</u>	<u>DAC</u> channel name. Default = spaces.
1346	sDACChannelUnits(0-3)	8 <u>char</u>	The user units for this <u>DAC</u> channel. Default = spaces.
1378	fDACScaleFactor(0-3)	<u>float</u>	<u>DAC</u> channel gain (user units / V at <u>DAC</u> ). Default = 1
1394	fDACHoldingLevel(0-3)	<u>float</u>	<u>DAC</u> channel holding level (user units). Default = 0.
1410	nSignalType	<u>short</u>	Type of <u>signal conditioner</u> that was used. 0 = None; 1 = CyberAmp 320/380.
1412	SUnused004	10 <u>char</u>	Unused.

## # \$ K + Synchronous timer outputs

(14 bytes).

Synchronous timer outputs were dropped from pCLAMP version 6. These parameters have been kept in the ABF 1.0 specification for compatibility when reading old data files. New programs should write zeros to all of the parameters in this group.

Offset	Header Entry Name	Type	Description
1422	nOUTEnable	<a href="#">short</a>	Enable synchronous timer outputs: 0 = No; 1 = Yes.
1424	nSampleNumberOUT1	<a href="#">short</a>	Sample number for pulse on synchronous timer OUT #1.
1426	nSampleNumberOUT2	<a href="#">short</a>	Sample number for pulse on synchronous timer OUT #2.
1428	nFirstEpisodeOUT	<a href="#">short</a>	First <a href="#">sweep</a> at which synchronous timer OUT #1 and #2 fire.
1430	nLastEpisodeOUT	<a href="#">short</a>	Last sweep at which synchronous timer OUT #1 and #2 fire.
1432	nPulseSamplesOUT1	<a href="#">short</a>	Duration and polarity of pulse on synchronous timer OUT #1 ( <a href="#">DAC</a> samples).
1434	nPulseSamplesOUT2	<a href="#">short</a>	Duration and polarity of pulse on synchronous timer OUT #2 ( <a href="#">DAC</a> samples).

---

# Synchronous\_timer\_outputs

\$ Synchronous timer outputs

K Synchronous timer  
outputs;nOUTEnable;nSampleNumberOUT1;nSampleNumberOUT2;nFirstEpisodeOUT;nLastEpisodeOUT;n  
PulseSamplesOUT1;nPulseSamplesOUT2;

+ ABFFUNC

## # \$ K + Epoch Waveform and Pulses

(184 bytes).

Description of epoch waveform properties. See Extended Epoch Waveform and Pulse for new parameters.

Offset	Header Entry Name	Type	Description
1436	nDigitalEnable	<a href="#">short</a>	Enable digital outputs: 0 = No; 1 = Yes.
1438	_nWaveformSource	<a href="#">short</a>	Analog waveform source: 0 = Disable; 1 = Generate waveform from epoch definitions; 2 = Generate waveform from a <a href="#">DAC</a> file.  <b>Discontinued:</b> Use <a href="#">nWaveformSource</a>
1440	nActiveDACChannel	<a href="#">short</a>	Active <a href="#">DAC</a> channel, i.e. the one used for waveform generation.  Redundant in ABF v1.8 except for control of the active digital output. To control each DAC waveform use <a href="#">nWaveformEnable</a>
1442	_nInterEpisodeLevel	<a href="#">short</a>	Inter- <a href="#">sweep</a> holding level: 0 = Use holding level; 1 = Use last epoch amplitude.  <b>Discontinued:</b> Use <a href="#">nInterEpisodeLevel</a>
1444	_nEpochType(0-9)	<a href="#">short</a>	Epoch type: 0 = Disabled; 1 = Step; 2 = Ramp.  <b>Discontinued:</b> Use <a href="#">nEpochType</a>
1464	_fEpochInitLevel(0-9)	<a href="#">float</a>	Epoch initial level (user units).  <b>Discontinued:</b> Use <a href="#">fEpochInitLevel</a>
1504	_fEpochLevelInc(0-9)	<a href="#">float</a>	Epoch level increment (user units).  <b>Discontinued:</b> Use <a href="#">fEpochLevelInc</a>
1544	_nEpochInitDuration(0-9)	<a href="#">short</a>	Epoch initial duration (in <a href="#">sequence</a> counts).  <b>Discontinued:</b> Use <a href="#">lEpochInitDuration</a>
1564	_nEpochDurationInc(0-9)	<a href="#">short</a>	Epoch duration increment (in sequence counts)  <b>Discontinued</b> <a href="#">lEpochDurationInc</a>
1584	nDigitalHolding	<a href="#">short</a>	Holding value for digital output.
1586	nDigitalInterEpisode	<a href="#">short</a>	Inter- <a href="#">sweep</a> digital holding value: 0 = Use holding value; 1 = Use last epoch value.
1588	nDigitalValue(0-9)	<a href="#">short</a>	Epoch value for digital output (0..9).
1610	sUnavailable1608	<a href="#">4char</a>	This parameter was used by CLAMPFIT 6.0, <b>do not reuse</b> . Was fWaveformOffset.  Offset (in active <a href="#">DAC</a> user units) in <a href="#">instrument</a> command pathway, usually due to a non-zero holding potential or current.
1612	nDigitalDACChannel	<a href="#">short</a>	Not used.
1614	sUnused005	<a href="#">6char</a>	Unused.

---

# Epoch\_Waveform\_and\_Pulses

\$ Epoch Waveform and Pulses

K Epoch Waveform and Pulses;nDigitalEnable;nWaveformSource;nActiveDACChannel;nInterEpisodeLevel;nEpochType;fEpochInitLevel;fEpochLevelInc;nEpochInitDuration;nEpochDurationInc;nDigitalHolding;nDigitalInterEpisode;nDigitalValue

+ ABFFUNC

## # \$ K + DAC Output File

(98 bytes).

Description of DAC File properties. See Extended [DAC](#) File for new parameters.

Offset	Header Entry Name	Type	Description
1620	_fDACFileScale	<a href="#">float</a>	Scaling factor to apply to <a href="#">DAC</a> file contents. <b>Discontinued:</b> Use <a href="#">fDACFileScale</a>
1624	_fDACFileOffset	<a href="#">float</a>	Offset (in user units) to apply to <a href="#">DAC</a> file contents. <b>Discontinued:</b> Use <a href="#">fDACFileOffset</a>
1628	sUnused006	<a href="#">2char</a>	Unused.
1630	_nDACFileEpisodeNum	<a href="#">short</a>	Sweep (or column) number to replay from waveform file: -1 = all except the first (which is skipped), repeating last if necessary; 0 = all sweeps; N = sweep number. <b>Discontinued:</b> Use <a href="#">IDACFileEpisodeNum</a>
1632	_nDACFileADCNum	<a href="#">short</a>	Logical <a href="#">ADC</a> channel number to replay from waveform file. <b>Discontinued:</b> Use <a href="#">nDACFileADCNum</a>
1634	_sDACFilePath	84 <a href="#">char</a>	Drive and directory for sDACFileName. <b>Discontinued:</b> Use <a href="#">sDACFilePath</a>

---

# DAC\_Output\_File

\$ DAC Output File

K DAC Output

File;fDACFileScale;fDACFileOffset;nDACFileEpisodeNum;nDACFileADCNum;sDACFileName;sDACFilePath  
;

+ ABFFUNC

## # \$ K + Conditioning pulse train

(44 bytes).

This section is replaced by Extended Conditioning Pulse Train variables. Also Conditioning Pulse Train has been renamed Pre-sweep Train in the application but variable names remain the same.

Offset	Header Entry Name	Type	Description
1718	_nConditEnable	<a href="#">short</a>	Conditioning pulse train activation status: 0 = Disable; 1 = Enable.  <b>Discontinued:</b> Use <a href="#">nConditEnable</a>
1720	_nConditChannel	<a href="#">short</a>	<a href="#">DAC</a> channel used for conditioning pulses.  <b>Discontinued.</b>
1722	_lConditNumPulses	<a href="#">long</a>	Number of pulses in conditioning pulse train.  <b>Discontinued:</b> Use <a href="#">lConditNumPulses</a>
1726	_fBaselineDuration	<a href="#">float</a>	A single pulse in the conditioning train consists of a baseline followed by a step. This parameter is the baseline duration in ms.  <b>Discontinued:</b> Use <a href="#">fBaselineDuration</a>
1730	_fBaselineLevel	<a href="#">float</a>	Baseline level (user units).  <b>Discontinued:</b> Use <a href="#">fBaselineLevel</a>
1734	_fStepDuration	<a href="#">float</a>	Step duration (ms).  <b>Discontinued:</b> Use <a href="#">fStepDuration</a>
1738	_fStepLevel	<a href="#">float</a>	Step level (user units).  <b>Discontinued:</b> Use <a href="#">fStepLevel</a>
1742	_fPostTrainPeriod	<a href="#">float</a>	At the end of the conditioning train there is a post-train steady-state output. This parameter is the post-train duration in ms.  <b>Discontinued:</b> Use <a href="#">fPostTrainPeriod</a>
1746	_fPostTrainLevel	<a href="#">float</a>	Post-train level (user units).  <b>Discontinued:</b> Use <a href="#">fPostTrainLevel</a>
1750	sUnused007	12 <a href="#">char</a>	Unused.

---

# Conditioning\_pulse\_train

\$ Conditioning pulse train

K Conditioning pulse train;nConditEnable;nConditChannel;lConditNumPulses;fBaselineDuration;fBaselineLevel;fStepDuration;fStepLevel;fPostTrainPeriod;fPostTrainLevel;

+ ABFFUNC

## # \$ K + Variable Parameter User List

(82 bytes).

This section is replaced by the Extended Variable Parameter User List.

Offset	Header Entry Name	Type	Description
1762	_nParamToVary	<a href="#">short</a>	Holds the index of the parameter that varies from <a href="#">sweep</a> to sweep in one <a href="#">run</a> .  <b>Discontinued:</b> Use <a href="#">nULParamToVary</a>
1764	_sParamValueList	80 <a href="#">char</a>	List of comma-separated values. If the number of entries in the list is fewer than the requested number of <a href="#">sweeps</a> , the last list value is re-used. If there are more values in the list, the excess list values are ignored.  <b>Discontinued:</b> Use <a href="#">sULParamValueList</a>

---

# Variable\_Parameter\_User\_List

\$ Variable Parameter User List

K Variable Parameter User List;nParamToVary;sParamValueList

+ ABFFUNC



## # \$ K + Statistics measurement

(36 bytes).

This section is replaced by the Extended Statistics Measurement.

Offset	Header Entry Name	Type	Description
1844	_nAutopeakEnable	<a href="#">short</a>	Autopeak activation status: 0 = Disabled; 1 = Enabled.  <b>Discontinued:</b> Use <a href="#">nStatsEnable</a>
1846	_nAutopeakPolarity	<a href="#">short</a>	-1 = search for negative peaks; 0 = search for absolute peak; 1 = search for positive peaks.  <b>Discontinued:</b> Use <a href="#">nStatsChannelPolarity</a>
1848	_nAutopeakADCNum	<a href="#">short</a>	Physical <a href="#">ADC</a> channel number to measure.  <b>Discontinued:</b> Use bit flag <a href="#">nStatsActiveChannels</a> for selecting multiple channels
1850	_nAutopeakSearchMode	<a href="#">short</a>	Search mode: 0-9 = epoch A-J; -1 = All; -2 = Use specified region.  <b>Discontinued:</b> Use <a href="#">nStatsSearchMode</a>
1852	_lAutopeakStart	<a href="#">long</a>	Start of specified statistics region. Only valid if nAutopeakSearchMode = -2.  <b>Discontinued:</b> Use <a href="#">lStatsStart</a>
1856	_AutopeakEnd	<a href="#">long</a>	End of specified statistics region. Only valid if nAutopeakSearchMode = -2.  <b>Discontinued:</b> Use <a href="#">lStatsEnd</a>
1860	_nAutopeakSmoothing	<a href="#">short</a>	Number of samples averaged in boxcar smoothing window when looking for a peak.  <b>Discontinued:</b> Use <a href="#">nStatsSmoothing</a>
1862	_nAutopeakBaseline	<a href="#">short</a>	Baseline for statistics measurements: -2 = None; -1 = Average the first (INumSamplesPerEpisode/64) samples; 0-9 = epoch A-J  <b>Discontinued:</b> Use <a href="#">nStatsBaseline</a>
1864	_nAutopeakAverage	<a href="#">short</a>	0 = search average and data; 1 = search average only. <b>Discontinued.</b>
1866	sUnavailable1866	<a href="#">2char</a>	Was nAutopeakSaveStrategy, use nStatisticsSaveStrategy
1868	_lAutopeakBaselineStart	<a href="#">long</a>	Start of baseline used in statistics measurements.  <b>Discontinued:</b> Use <a href="#">lStatsBaselineStart</a>
1872	_lAutopeakBaselineEnd	<a href="#">long</a>	End of baseline used in statistics measurements.  <b>Discontinued:</b> Use <a href="#">lStatsBaselineEnd</a>
1876	_lAutopeakMeasurements	<a href="#">long</a>	Bit mask indicating which statistics measurements to perform.  <b>Discontinued:</b> Use <a href="#">lStatsMeasurements</a>

---

# Autopeak\_measurement

\$ Autopeak measurement

K Autopeak

measurement;nAutopeakEnable;nAutopeakPolarity;nAutopeakADCNum;nAutopeakSearchMode;lAutopeakStart;lAutopeakEnd;nAutopeakSmoothing;nAutopeakBaseline;nAutopeakAverage;

+ ABFFUNC

## # \$ K + Channel Arithmetic

(52 bytes).

Several standard arithmetic expressions are supported:

Expression #0 (general purpose)

Result = (K1\*A + K2) <op> (K3\*B + K4)

Expression #1 (ratio dyes)

Result = (K1\*R + K2) <op> (K3\*R + K4)

where

$R = (A + K5)/(B + K6)$

In both cases:

A and B are [ADC](#) channel numbers. A = B is allowed.

K1 .. K6 are constants.

<op> is an arithmetic operator.

Offset	Header Entry Name	Type	Description
1880	nArithmeticEnable	<a href="#">short</a>	Arithmetic activation status: 0 = Disable; 1 = Enable
1882	fArithmeticUpperLimit	<a href="#">float</a>	Display upper limit for arithmetic Results channel.
1886	fArithmeticLowerLimit	<a href="#">float</a>	Display lower limit for arithmetic Results channel.
1890	nArithmeticADCNumA	<a href="#">short</a>	Physical ADC channel number used for A.
1892	nArithmeticADCNumB	<a href="#">short</a>	Physical ADC channel number used for B.
1894	fArithmeticK1	<a href="#">float</a>	Numeric constant K1
1898	fArithmeticK2	<a href="#">float</a>	Numeric constant K2
1902	fArithmeticK3	<a href="#">float</a>	Numeric constant K3
1906	fArithmeticK4	<a href="#">float</a>	Numeric constant K4
1910	sArithmeticOperator	<a href="#">2char</a>	Arithmetic operator: "+", "-", "*", or "/"
1912	sArithmeticUnits	<a href="#">8char</a>	Units for Results channel.
1920	fArithmeticK5	<a href="#">float</a>	Numeric constant K5
1924	fArithmeticK6	<a href="#">float</a>	Numeric constant K6
1928	NArithmeticExpression	<a href="#">short</a>	The expression to evaluate: N = Expression #N.
1930	sUnused1930	<a href="#">2char</a>	Unused.

---

# Channel\_Arithmetic

\$ Channel Arithmetic

<sup>K</sup> Channel

Arithmetic;nArithmeticEnable;fArithmeticUpperLimit;fArithmeticLowerLimit;nArithmeticADCNumA;nArithmeticADCNumB;fArithmeticK1;fArithmeticK2;fArithmeticK3;fArithmeticK4;sArithmeticOperator;sArithmeticUnits;fArithmeticK5;fArithmeticK6;nArithmeticExpression;

+ ABFFUNC

## # \$ K + On-line Subtraction

(34 bytes).

Parts of this section have been replaced by the Extended On-line subtraction.

Offset	Header Entry Name	Type	Description
1932	_nPNEnable	<a href="#">short</a>	P/N activation status: 0 = Disable; 1 = Enable  <b>Discontinued:</b> Use <a href="#">nPNEnable</a>
1934	nPNPosition	<a href="#">short</a>	P/N subpulse execution: 0 = Before waveform; 1 = After waveform.
1936	_nPNPolarity	<a href="#">short</a>	-1 = opposite polarity to the waveform; 1 = same polarity as the waveform.  <b>Discontinued:</b> Use <a href="#">nPNPolarity</a>
1938	nPNNumPulses	<a href="#">short</a>	Number of P/N subpulses.
1940	_nPNADCNum	<a href="#">short</a>	Physical <a href="#">ADC</a> channel number for P/N subtraction.  <b>Discontinued:</b> Use <a href="#">nPNADCNum</a>
1942	_fPNHoldingLevel	<a href="#">float</a>	Subpulse holding level.  <b>Discontinued:</b> Use <a href="#">fPNHoldingLevel</a>
1946	fPNSettlingTime	<a href="#">float</a>	Settling time at subpulse holding level (ms). Wait this interval after changing to subpulse holding level before starting subpulses; wait this interval after returning to normal holding level before starting epochs.
1950	fPNInterpulse	<a href="#">float</a>	Timing interval between the start of one subpulse and the start of the next (ms).
1954	sUnused009	12 <a href="#">char</a>	Unused.

---

# On\_line\_Subtraction

\$ On-line Subtraction

K On-line Subtraction

nPNEnable;nPNPosition;nPNPolarity;nPNNumPulses;nPNADCNum;fPNHoldingLevel;fPNSettlingTime;fPNInterpulse;

+ ABFFUNC

## # \$ K + Miscellaneous Parameters

(82 bytes).

Offset	Header Entry Name	Type	Description
1966	_nListEnable	<a href="#">short</a>	Parameter list activation status: 0 = Disable; 1 = Enable.  <b>Discontinued:</b> Use <a href="#">nULEnable</a>
1968	nBellEnable(0-1)	<a href="#">short</a>	Auditory tone activation status: 0 = Disable; 1 = Enable
1972	nBellLocation(0-1)	<a href="#">short</a>	Location of bell relative to <a href="#">trial</a> : 0 = Before; 1 = After
1976	nBellRepetitions(0-1)	<a href="#">short</a>	Number of sounds to produce.
1980	nLevelHysteresis	<a href="#">short</a>	Amount of level hysteresis to use when detecting events. This is the amount that the data has to go past the trigger level before it is considered triggered. Re-arming of the trigger level is always done at the actual nominated trigger level. (See fTriggerThreshold)
1982	lTimeHysteresis	<a href="#">long</a>	Amount of time hysteresis to use when detecting events. This is the number of samples that have to be blow the trigger point before the trigger is said to be rearmed.
1986	nAllowExternalTags	<a href="#">short</a>	0 = Do not scan for external tags during acquisition. 1 = Scan for external tags.
1988	nLowpassFilterType(0-15)	<a href="#">char</a>	Type of Low Pass filter for each <a href="#">ADC</a> channel: 0 = None; 1 = External; 2 = Simple RC; 3 = Bessell; 4 = Butterworth.
2004	nHighpassFilterType(0-15)	<a href="#">char</a>	Type of High Pass filter for each ADC channel: 0 = None; 1 = External; 2 = Simple RC; 3 = Bessell; 4 = Butterworth.
2020	nAverageAlgorithm	<a href="#">short</a>	Algorithm used for calculating averages: 0 = Cumulative Averaging; 1 = Most Recent Averaging (uses fAverageWeighting below)
2022	fAverageWeighting	<a href="#">float</a>	Weighting Factor for Most Recent Averaging. This is the proportion of the incoming <a href="#">sweep</a> to include in the average.
2026	nUndoPromptStrategy	<a href="#">short</a>	Strategy for Prompting to create an Undo file: 0 = On Abort; 1 = Always;
2028	nTrialTriggerSource	<a href="#">short</a>	Trigger source for start <a href="#">trial</a> : -3 = Spacebar; -2 = External Trigger; -1 = None
2030	nStatisticsDisplayStrategy	<a href="#">short</a>	Strategy for displaying statistics: 0 = Display Statistics; 1 = Do Not display Statistics
2032	nExternalTagType	<a href="#">short</a>	Type of External Tag: 0 = Time Tag; 1 = External Tag; 2 = External Tag (from BNC input) 3 = Voice Tag (from audio input); 4 = New File Tag
2034	lHeaderSize	<a href="#">long</a>	Total size of the header. Currently 6144 bytes.

---

# Miscellaneous\_parameters

\$ Miscellaneous parameters

K Miscellaneous parameters

nListEnable;nBellEnable;nBellLocation;nBellRepetitions;nLevelHysteresis;lTimeHysteresis;nAllowExternalTags;

+ ABFFUNC

2038	dFileDuration	<u>double</u>	Not used.
2046	nStatisticsClearStrategy	<u>short</u>	Strategy for clearing statistics from display: 0 = Do not clear statistics; 1 = Clear Statistics.

# \$ K + **Extended File Structure**

(26 bytes).

Offset	Header Entry Name	Type	Description
2048	IDACFilePtr(0-1)	<a href="#">long</a>	Pointer to location of <a href="#">DAC</a> stimulus file data in physical <a href="#">DAC</a> channel order.
2056	IDACFileNumEpisodes(0-1)	<a href="#">long</a>	Number of episodes for each <a href="#">DAC</a> stimulus file in physical <a href="#">DAC</a> channel order.
2064	sUnused010	10 <a href="#">char</a>	Unused.

---

# Extended\_File\_Structure

\$ Extended File Structure

K Extended File Structure

IDACFilePtr;IDACFileNumEpisodes;

+ ABFFUNC

## # \$ K + Extended Multi-Channel Information

(62 bytes)

Offset	Header Entry Name	Type	Description
2074	fDACCcalibrationFactor(0-1)	<a href="#">float</a>	Calibration factor for the current digitizer in physical <a href="#">DAC</a> channel order.
2090	fDACCcalibrationOffset(0-1)	<a href="#">float</a>	Calibration offset for the current digitizer in physical <a href="#">DAC</a> channel order.
2106	sUnused011	30 <a href="#">char</a>	Unused.

---

# Extended\_MultiChannel\_Information

\$ Extended MultiChannel Information

K Extended MultiChannel Information

fDACCcalibrationFactor;fDACCcalibrationOffset;

+ ABFFUNC

## # \$ K + Train Parameters

(160 bytes).

Offset	Header Entry Name	Type	Description
2136	IEpochPulsePeriod(0-1)(0-9)	<a href="#">long</a>	Train period in physical <a href="#">DAC</a> channel order then epoch order (samples).
2216	IEpochPulseWidth(0-1)(0-9)	<a href="#">long</a>	Train pulse width in physical <a href="#">DAC</a> channel order then epoch order (samples).

---

# Train\_Parameters

\$ Train Parameters

<sup>K</sup> Train Parameters

IEpochPulsePeriod;IEpochPulseWidth;

<sup>+</sup> ABFFUNC



## # \$ K + Extended Epoch Waveform and Pulses

(412 bytes).

Offset	Header Entry Name	Type	Description
2296	nWaveformEnable(0-1)	<a href="#">short</a>	Analog waveform output status for each physical <a href="#">DAC</a> channel: 0 = Disabled, 1 = Enabled.
2300	nWaveformSource(0-1)	<a href="#">short</a>	Analog waveform source in physical <a href="#">DAC</a> channel order: 0 = Disable; 1 = Generate waveform from epoch definitions; 2 = Generate waveform from a <a href="#">DAC</a> file.
2304	nInterEpisodeLevel(0-1)	<a href="#">short</a>	Inter-sweep holding level in physical <a href="#">DAC</a> channel order: 0 = Use holding level; 1 = Use last epoch amplitude.
2308	nEpochType(0-1)(0-9)	<a href="#">short</a>	Epoch type in physical <a href="#">DAC</a> channel order then epoch order: 0 = Disabled; 1 = Step; 2 = Ramp.
2348	fEpochInitLevel(0-1)(0-9)	<a href="#">float</a>	Epoch initial level in physical <a href="#">DAC</a> channel order then epoch order (user units).
2428	fEpochLevelInc(0-1)(0-9)	<a href="#">float</a>	Epoch level increment in physical <a href="#">DAC</a> channel order then epoch order (user units).
2508	lEpochInitDuration(0-1)(0-9)	<a href="#">long</a>	Epoch initial duration in physical <a href="#">DAC</a> channel order then epoch order (in <a href="#">sequence</a> counts).
2588	lEpochDurationInc(0-1)(0-9)	<a href="#">long</a>	Epoch duration increment in physical <a href="#">DAC</a> channel order then epoch order (in <a href="#">sequence</a> counts)
2668	nDigitalTrainValue(0-9)	<a href="#">short</a>	Epoch value for digital train output in epoch order. 0000 = Disabled; 0*000 = Generates digital train on bit 3. Train period and pulse width can be controlled by the user list.  <div style="text-align: right;">           ASCII digital train pattern 0000 :0x00000000             000* :0x00000001             00*0 :0x00000002             0*00 :0x00000004             *000 :0x00000008         </div>
2688	nDigitalTrainActiveLogic	<a href="#">short</a>	Digital train polarity. 0 = train starts with LOW pulse. 1 = trains starts with HIGH pulse.
2690	sUnused012	18 <a href="#">char</a>	Unused.

---

# Extended\_Epoch\_Waveform\_and\_Pulses

\$ Extended Epoch Waveform and Pulses

K Extended Epoch Waveform and Pulses

nWaveformEnable;nWaveformSource;nInterEpisodeLevel;nEpochType;fEpochInitLevel;fEpochLevelInc;lEpochInitDuration;lEpochDurationInc;nDigitalTrainValue;nDigitalTrainActiveLogic;

+ ABFFUNC

## # \$ K + Extended DAC Output File

(552 bytes).

Offset	Header Entry Name	Type	Description
2708	fDACFileScale(0-1)	<a href="#">float</a>	Scaling factor to apply to <a href="#">DAC</a> file contents in physical <a href="#">DAC</a> channel order.
2716	fDACFileOffset(0-1)	<a href="#">float</a>	Offset (in user units) to apply to <a href="#">DAC</a> file contents in physical <a href="#">DAC</a> channel order..
2724	IDACFileEpisodeNum(0-1)	<a href="#">long</a>	Sweep (or column) number to replay from waveform file in physical <a href="#">DAC</a> channel order: -1 = all except the first (which is skipped), repeating last if necessary; 0 = all sweeps; N = sweep number.
2732	nDACFileADCNum(0-1)	<a href="#">short</a>	Logical <a href="#">ADC</a> channel number to replay from waveform file in physical <a href="#">DAC</a> channel order.
2736	sDACFilePath(0-1)	256 <a href="#">char</a>	Drive and directory for sDACFileName in physical <a href="#">DAC</a> channel order.
3248	sUnused013	12 <a href="#">char</a>	Unused.

---

# Extended\_DAC\_Output File

\$ Extended DAC Output File

K Extended DAC Output File

fDACFileScale;fDACFileOffset;IDACFileEpisodeNum;nDACFileADCNum;sDACFilePath;

+ ABFFUNC

## # \$ K + Extended Pre-sweep (Conditioning) Pulse Train

(100 bytes).

Offset	Header Entry Name	Type	Description
3260	nConditEnable(0-1)	<a href="#">short</a>	Pre-sweep pulse train activation status in physical <a href="#">DAC</a> channel order: 0 = Disable; 1 = Enable.
3264	lConditNumPulses(0-1)	<a href="#">long</a>	Number of pre-sweep train pulses in physical <a href="#">DAC</a> channel order.
3272	fBaselineDuration(0-1)	<a href="#">float</a>	A single pulse in the pre-sweep train consists of a baseline followed by a step. This parameter is the baseline duration in ms in physical <a href="#">DAC</a> channel order.
3280	fBaselineLevel(0-1)	<a href="#">float</a>	Baseline level in physical <a href="#">DAC</a> channel order (user units).
3288	fStepDuration(0-1)	<a href="#">float</a>	Step duration in physical <a href="#">DAC</a> channel order (ms).
3296	fStepLevel(0-1)	<a href="#">float</a>	Step level in physical <a href="#">DAC</a> channel order (user units).
3304	fPostTrainPeriod(0-1)	<a href="#">float</a>	At the end of the conditioning train there is a post-train steady-state output. This parameter is the post-train duration in ms in physical <a href="#">DAC</a> channel order.
3312	fPostTrainLevel(0-1)	<a href="#">float</a>	Post-train level in physical <a href="#">DAC</a> channel order (user units).
3320	sUnused014	40 <a href="#">char</a>	Unused.

---

# Extended\_Presweep\_(Conditioning)\_Pulse\_Train

\$ Extended Presweep (Conditioning) Pulse Train

K Extended Presweep (Conditioning) Pulse Train;

fDACFileScale;fDACFileOffset;lDACFileEpisodeNum;nDACFileADCNum;sDACFilePath;

+ ABFFUNC

## # \$ K + Extended Variable Parameter User List

(1096 bytes).

Offset	Header Entry Name	Type	Description
3360	nULEnable(0-3)	<a href="#">short</a>	Parameter list activation status in physical <a href="#">DAC</a> channel order: 0 = Disable; 1 = Enable.
3368	nULParamToVary(0-3)	<a href="#">short</a>	Holds the index of the parameter that varies from <a href="#">sweep</a> to sweep in one <a href="#">run</a> in physical <a href="#">DAC</a> channel order.
			CONDITNUMPULSES 0
			CONDITBASELINEDURATION 1
			CONDITBASELINELEVEL 2
			CONDITSTEPDURATION 3
			CONDITSTEPLEVEL 4
			CONDITPOSTTRAINDURATION 5
			CONDITPOSTTRAINLEVEL 6
			EPISODESTARTTOSTART 7
			INACTIVEHOLDING 8
			DIGITALINTEREPIISODE 9
			PNNUMPULSES 10
			PARALLELVALUE(0-9) 11-20
			EPOCHINITLEVEL(0-9) 21-30
			EPOCHINITDURATION(0-9) 31-40
			EPOCHTRAINPERIOD(0-9) 41-50
			EPOCHTRAINPULSEWIDTH(0-9) 51-60
3376	sULParamValueList(0-3)	256 <a href="#">char</a>	List of comma-separated values in physical <a href="#">DAC</a> channel order. If the number of entries in the list is fewer than the requested number of <a href="#">sweeps</a> , the last list value is re-used but only if nULRepeat = 0. If there are more values in the list, the excess list values are ignored.
4400	nULRepeat(0-3)	<a href="#">short</a>	Repeat the list when the current sweep exceeds the number of entries in the list. 0 = Disable, 1 = Repeat the list.
4408	sUnused015	48 <a href="#">char</a>	Unused.

---

# Extended\_Variable\_Parameter\_User\_List

\$ Extended Variable Parameter User List

K Extended Variable Parameter User List

nULEnable;nULParamToVary;sULParamValueList;nULRepeat;

+ ABFFUNC

## # \$ K + Extended On-Line Subtraction

(56 bytes).

Offset	Header Entry Name	Type	Description
4456	nPNEnable(0-1)	<a href="#">short</a>	P/N activation status in physical <a href="#">DAC</a> channel order: 0 = Disable; 1 = Enable
4460	nPNPolarity(0-1)	<a href="#">short</a>	Polarity of P/N waveform in physical <a href="#">DAC</a> channel order: -1 = opposite polarity to the waveform; 1 = same polarity as the waveform.
4464	nPNADCNum(0-1)	<a href="#">short</a>	Physical <a href="#">ADC</a> channel number for P/N subtraction in physical <a href="#">DAC</a> channel order.
4468	fPNHoldingLevel(0-1)	<a href="#">float</a>	Subpulse holding level in physical <a href="#">DAC</a> channel order.
4476	sUnused016	36 <a href="#">char</a>	Unused.

---

# Extended\_On\_Line\_Subtraction

\$ Extended On-Line Subtraction

K Extended On-Line Subtraction

nPNEnable;nPNPolarity;nPNADCNum;fPNHoldingLevel;

+ ABFFUNC

## # \$ K + Extended Environmental Information

(898 bytes).

Offset	Header Entry Name	Type	Description
4512	nTelegraphEnable(0-15)	<a href="#">short</a>	Enable storage of telegraph information in physical <a href="#">ADC</a> channel order: 0 = Disabled; 1 = Automatic; 2 = Manual.
4544	nTelegraphInstrument(0-15)	<a href="#">short</a>	Telegraphed <a href="#">instrument</a> in physical <a href="#">ADC</a> channel order. (IHS-1 telegraphs are not supported.) Note, for most programs this is an information-only field. For example, in Clampex the telegraphed instrument is chosen as a configuration item and copied into this field.
			Unknown instrument 0
			Axopatch-1 with CV-4-1/100 1
			Axopatch-1 with CV-4-0.1/100 2
			Axopatch-1B(inv.) CV-4-1/100 3
			Axopatch-1B(inv) CV-4-0.1/100 4
			Axopatch 200 with CV 201 5
			Axopatch 200 with CV 202 6
			GeneClamp 7
			Dagan 3900 8
			Dagan 3900A 9
			Dagan CA-1 Im=0.1 10
			Dagan CA-1 Im=1.0 11
			Dagan CA-1 Im=10 12
			Warner OC-725 13
			Warner OC-725C 14
			AxoPatch 200B 15
			Dagan PC-ONE Im=0.1 16
			Dagan PC-ONE Im=1.0 17
			Dagan PC-ONE Im=10 18
			Dagan PC-ONE Im=100 19
			Warner BC-525C 20
			Warner PC-505 21
			Warner PC-501 22
			Dagan CA-1 Im=0.05 23

---

# Extended\_Environmental\_Information

\$ Extended Environmental Information

K Extended Environmental Information

nTelegraphEnable;nTelegraphInstrument;fTelegraphAdditGain;fTelegraphFilter;fTelegraphMembraneCap;nTelegraphMode;nManualTelegraphStrategy;nAutoAnalyseEnable;sAutoAnalysisMacroName;sProtocolPath;sFileComment;

+ ABFFUNC

			MultiClamp 700	24	
			Turbo Tec	25	
			OpusXpress 6000	26	
4576	fTelegraphAdditGain(0-15)	<u>float</u>	Additional gain multiplier of telegraphed <u>instrument</u> in physical <u>ADC</u> channel order. (Optionally autosampled by some acquisition programs.) (Default = 1)		
4640	fTelegraphFilter(0-15)	<u>float</u>	Lowpass filter cutoff frequency of telegraphed <u>instrument</u> in <u>ADC</u> physical channel order. (Optionally autosampled by some acquisition programs.) (Default = 100000)		
4704	fTelegraphMembraneCap(0-15)	<u>float</u>	Patch-clamp membrane capacitance compensation of telegraphed <u>instrument</u> in physical <u>ADC</u> channel order. (Optionally autosampled by some acquisition programs.)		
4768	nTelegraphMode(0-15)	<u>short</u>	Telegraph mode in physical <u>ADC</u> channel order: 1 = Voltage Clamp; 2 = Current Clamp; 4 = Current Clamp Zero (Automatically set by MultiClamp <u>instrument</u> )		
4800	nTelegraphDACScaleFactor Enable(0-3)	<u>short</u>	Determines whether fDACScaleFactor was telegraphed: 1 = telegraphed; 0 = not telegraphed		
4808	sUnused016a	24 <u>char</u>	Unused.		
4832	nAutoAnalyseEnable	<u>short</u>	Discontinued.		
4834	sAutoAnalysisMacroName	64 <u>char</u>	Discontinued.		
4898	sProtocolPath	256 <u>char</u>	The name and path of the protocol		
5154	sFileComment	128 <u>char</u>	128 byte ASCII comment string.		
5282	sUnused017	128 <u>char</u>	Unused.		

## # \$ K + Extended Statistics Measurements

(388 bytes).

Offset	Header Entry Name	Type	Description
5410	nStatsEnable	<u>short</u>	Statistics activation status: 0 = Disabled; 1 = Enabled.
5412	nStatsActiveChannels	<u>unsigned short</u>	Bit flag for selecting statistics in each physical <u>ADC</u> channel.
			Channel 0: 0x0001
			Channel 1: 0x0002
			Channel 2: 0x0004
			Channel 3: 0x0008
			Channel 4: 0x0010
			Channel 5: 0x0020
			Channel 6: 0x0040
			Channel 7: 0x0080
			Channel 8: 0x0100
			Channel 9: 0x0200
			Channel 10: 0x0400
			Channel 11: 0x0800
			Channel 12: 0x1000
			Channel 13: 0x2000
			Channel 14: 0x4000
			Channel 15: 0x8000
			All channels: 0xFFFF
5414	nStatsSearchRegionFlags	<u>unsigned short</u>	Bit flag for selecting statistics in each region.
			Region 0: 0x0001
			Region 1: 0x0002
			Region 2: 0x0004
			Region 3: 0x0008
			Region 4: 0x0010
			Region 5: 0x0020
			Region 6: 0x0040
			Region 7: 0x0080

---

# Extended\_Statistics\_Measurements

\$ Extended Statistics Measurements

K Extended Statistics Measurements

nStatsEnable;nStatsActiveChannels;nStatsSearchRegionFlags;nStatsSelectedRegion;\_nStatsSearchMode;nStatsSmoothing;nStatsSmoothingEnable;nStatsBaseline;lStatsBaselineStart;lStatsBaselineEnd;lStatsMeasurements;lStatsStart;lStatsEnd;nRiseBottomPercentile;RiseTopPercentile;nDecayBottomPercentile;nDecayTopPercentile;nStatsChannelPolarity;nStatsSearchMode;

+ ABFFUNC



			All regions: 0x00FF
5416	nStatsSelectedRegion	<u>short</u>	The currently active statistics region: 0 = region 0; 1 = region 1 ... 7 = region 7.
5418	_nStatsSearchMode	<u>short</u>	Not used
5420	nStatsSmoothing	<u>short</u>	Number of samples averaged in boxcar smoothing window when looking for a peak.
5422	nStatsSmoothingEnable	<u>short</u>	Not used
5422	nStatsBaseline	<u>short</u>	Baseline for statistics measurements: -2 = None; -1 = Average the first (INumSamplesPerEpisode/64) samples; 0-9 = epoch A-J
5426	IStatsBaselineStart	<u>long</u>	Start of baseline used in statistics measurements (samples).
5430	IStatsBaselineEnd	<u>long</u>	End of baseline used in statistics measurements (samples)
5434	IStatsMeasurements(0-7)	<u>long</u>	Bit mask indicating which statistics measurements to perform in region order.

Peak: 0x00000001

Peak time: 0x00000002

Antipeak: 0x00000004

Antipeak Time: 0x00000008

Mean: 0x00000010

Standard Deviation: 0x00000020

Integral: 0x00000040

Max Rise Slope: 0x00000080

Max Rise Slope Time: 0x00000100

Max Decay Slope: 0x00000200

Max Decay Slope Time: 0x00000400

Rise Time: 0x00000800

Decay Time: 0x00001000

Half Width: 0x00002000

Baseline: 0x00004000

Rise Slope: 0x00008000

Decay Slope: 0x00010000

Region Slope: 0x00020000

All Measurements: 0x0002FFFF

5466	IStatsStart(0-7)	<u>long</u>	Start of specified statistics region in region order (Samples). Only valid if corresponding value in nStatsSearchMode = -2.
5498	IStatsEnd(0-7)	<u>long</u>	End of specified statistics region in region order (Samples). Only valid if corresponding value in nStatsSearchMode = -2.
5530	nRiseBottomPercentile(0-7)	<u>short</u>	Nominated percentage of the peak amplitude (relative to the baseline) that determines the left most point of region for measuring rise slope parameters: 0 = at baseline; 100 = at peak.

5546	nRiseTopPercentile(0-7)	<u>short</u>	Nominated percentage of the peak amplitude (relative to the baseline) that determines the right most point of region for measuring rise slope parameters. 0 = at baseline; 100 = at peak.
5562	nDecayBottomPercentile(0-7)	<u>short</u>	Nominated percentage of the peak amplitude (relative to the baseline) that determines the right most point of region for measuring decay slope parameters. 0 = at baseline; 100 = at peak.
5578	nDecayTopPercentile(0-7)	<u>short</u>	Nominated percentage of the peak amplitude (relative to the baseline) that determines the left most point of region for measuring decay slope parameters. 0 = at baseline; 100 = at peak.
5594	nStatsChannelPolarity(0-15)	<u>short</u>	Polarity of peaks to find in physical <u>ADC</u> channel order. -1 = search for negative peaks; 0 = search for absolute peak; 1 = search for positive peaks.
5626	nStatsSearchMode(0-7)	<u>short</u>	Search mode in region order: 0-9 = epoch A-J; -1 = All; -2 = Use specified region.
5642	sUnused018	156 <u>char</u>	Unused.

## # \$ K + Application Version Data

(16 bytes).

Offset	Header Entry Name	Type	Description
5798	nMajorVersion	<a href="#">short</a>	Major version of application: x.0.0.0
5800	nMinorVersion	<a href="#">short</a>	Minor version of application: 0.x.0.0
5802	nBugfixVersion	<a href="#">short</a>	Bug fix version of application: 0.0.x.0
5804	nBuildVersion	<a href="#">short</a>	Build version of application: 0.0.0.x
5806	sUnused019	<a href="#">8char</a>	Unused.

---

# Application\_Version\_Data

\$ Application Version Data

K Application Version Data

nMajorVersion;nMinorVersion;nBugfixVersion;nBuildVersion;

+ ABFFUNC

## # \$ K + LTP Protocol

(14 bytes).

Offset	Header Entry Name	Type	Description
5814	nLTPTType	<u>char</u>	Type of LTP protocol: 0 = None; 1 = Baseline; 2 = Induction/Conditioning
5816	nLTPUsageOfDAC(0-1)	<u>short</u>	Indicates how specified <u>DAC</u> is used: 0 = Presynaptic stimulus; Postsynaptic stimulus.
5820	nLTPPresynapticPulses(0-1)	<u>short</u>	The number pulses in presynaptic pulse train.
5824	sUnused020	<u>4char</u>	Unused.

---

# LTP\_Protocol

\$ LTP Protocol

K LTP Protocol

nLTPTType;nLTPUsageOfDAC;nLTPPresynapticPulses;

+ ABFFUNC

## # \$ K + Digidata 132x Trigger Out

(8 bytes).

Offset	Header Entry Name	Type	Description
5828	nDD132xTriggerOut	<a href="#">short</a>	The trigger enable for DD132x series digitizers
5830	sUnused021	<a href="#">6char</a>	Unused.

---

# Digidata\_132x\_Trigger\_Out

\$ Digidata 132x Trigger Out

K Digidata 132x Trigger Out

nDD132xTriggerOut;

+ ABFFUNC

## # \$ K + Epoch Resistance

(40 bytes).

Offset	Header Entry Name	Type	Description
5836	sEpochResistanceSignalName(0-1)	10 <u>char</u>	Name of the signal on which to measure epoch resistance in physical <u>DAC</u> channel order.
5856	nEpochResistanceState(0-1)	<u>short</u>	Epoch resistance status in physical <u>DAC</u> channel order: 0 = Disabled; 1 = Enabled.
5860	sUnused022	16 <u>char</u>	Unused.

---

# Epoch\_Resistance

\$ Epoch Resistance

K Epoch Resistance

sEpochResistanceSignalName;nEpochResistanceState;

+ ABFFUNC

## # \$ K + Alternating Episodic Mode

(58 bytes).

Offset	Header Entry Name	Type	Description
5876	nAlternateDACOutputState	<a href="#">short</a>	Alternating <a href="#">DAC</a> output status: 0 = Disabled; 1 = Enabled. For the first <a href="#">sweep</a> <a href="#">DAC</a> 0 is set by waveform 0 and <a href="#">DAC</a> 1 is set to holding level. For the second <a href="#">sweep</a> <a href="#">DAC</a> 1 is set by waveform 1 and <a href="#">DAC</a> 0 is set to holding level etc.
5878	nAlternateDigitalValue(0-9)	<a href="#">short</a>	Digital value held by the alternate or non-active digital output: The active digital output is set by <a href="#">nActiveDACChannel</a> .
5898	nAlternateDigitalTrainValue(0-9)	<a href="#">short</a>	Digital train value held by the alternate or non-active digital output: The active digital output is set by <a href="#">nActiveDACChannel</a> .
5918	nAlternateDigitalOutputState	<a href="#">short</a>	Alternating digital output status: 0 = Disabled; 1 = Enabled. . For the first <a href="#">sweep</a> the digital pattern is set by waveform 0. For the second <a href="#">sweep</a> the digital pattern is set by waveform 1 etc.
5920	sUnused023	14 <a href="#">char</a>	Unused.

---

# Alternating\_Episodic\_Mode

\$ Alternating Episodic Mode

K Alternating Episodic Mode

nAlternateDACOutputState;nAlternateDigitalValue;nAlternateDigitalTrainValue;nAlternateDigitalOutputState;

+ ABFFUNC

## # \$ K + Post-Processing Actions

(210 bytes).

Offset	Header Entry Name	Type	Description
5934	fPostProcessLowpassFilter(0-15)	<a href="#">float</a>	Post processing lowpass filter cutoff for single io channel data (Hz)
5998	nPostProcessLowpassFilterType(0-15)	<a href="#">char</a>	Post processing lowpass filter type
			None 0
			Combination 1
			Bessel 2
			Boxcar 3
			Butterworth 4
			Chebyshev 5
			Gaussian 6
			Rc 7
			Rc8 8
			Adaptive 9
6014	sUnused2048	130 <a href="#">char</a>	Unused.

TOTAL 6144 bytes

---

# Post\_Processing\_Actions

\$ Post-Processing Actions

K Post-Processing Actions

fPostProcessLowpassFilter;nPostProcessLowpassFilterType;

+ ABFFUNC



## # \$ K + The ABF Scope Config Section

If present, the ABF Scope Config section will contain one or more ABFScopeConfig structures describing the attributes of the scope windows used for data display during the data acquisition.

### The ABFScopeConfig Structure

Each ABFScopeConfig structure contains configuration information that describes the setup of a scope window. This structure in turn contains a structure (ABFLogFont) to define the font properties used to draw textual items such as tic labels, and an array of structures (of type ABFSignal), one for each [ADC](#) channel being acquired.

Offset	Header Entry Name	Type	Description
0	dwFlags	<a href="#">DWORD</a>	Flags that are meaningful to the scope.
4	rgbColor[ABF_SCOPECOLORS]	DWORD	Colors for the components of the scope.
44	fDisplayStart	<a href="#">float</a>	Start of the display area in ms.
48	fDisplayEnd	float	End of the display area in ms.
52	wScopeMode	<a href="#">WORD</a>	The display mode: 0=sweeps; 1=continuous.
54	bMaximized	char	TRUE = Scope parent is maximized.
55	bMinimized	char	TRUE = Scope parent is minimized.
56	xLeft	<a href="#">short</a>	Coordinate of the left edge.
58	yTop	<a href="#">short</a>	Coordinate of the top edge.
60	xRight	<a href="#">short</a>	Coordinate of the right edge.
62	yBottom	<a href="#">short</a>	Coordinate of the bottom edge.
64	LogFont	ABFLogFont	Description of current font.
104	TraceList[ABF_ADCCOUNT]	ABFSignal	List of traces in current use - see nADCNumChannels in the ABFFileHeader for the number of channels actually used.
648	nYAxisWidth	<a href="#">short</a>	Width of the YAxis region.
650	nTraceCount	<a href="#">short</a>	Number of traces described in TraceList. (this should always match nADCNumChannels in the ABFFileHeader)
652	nEraseStrategy	<a href="#">short</a>	Erase strategy: 0=Erase before each sweep; 1=Erase before each run; 2=Erase before each trial; 3=Do not erase
654	nDockState	<a href="#">short</a>	Docking state: 0=Not docked; 1=Top; 2=Left; 3=Right; 4=Bottom.

TOTAL 656 bytes

### The ABFLogFont Structure

The ABFLogFont structure is a subset of the Windows LogFont structure, containing information to describe the characteristics of the font to be used in the Scope window.

---

# The\_ABF\_Scope\_Config\_Section

\$ The ABF Scope Config Section

K The ABF Scope Config Section;ABFLogFont;ABFSignal;ABFScopeConfig

+ ABFFUNC

Offset	Header Entry Name	Type	Description
0	nHeight	<a href="#">short</a>	N.B. Height of the font in *points*.
2	nWeight	<a href="#">short</a>	MSWindows font weight value.
4	cPitchAndFamily	char	MSWindows pitch and family mask.
5	Unused[3]	3char	Unused space to maintain 4-byte packing.
8	szFaceName[ABF_FACESIZE]	char	Face name of the font.
<u>TOTAL 40 bytes</u>			

### The ABFSignal Structure

The ABFSignal structure describes the characteristics of a single data trace on the screen, corresponding to a particular acquired signal.

Offset	Header Entry Name	Type	Description
0	szName[ABF_ADCNAMELEN+2]	<a href="#">char</a>	ABF name length + '\0' + 1 for alignment.
12	nMxOffset	<a href="#">short</a>	Offset of the signal in the sampling sequence.
14	RgbColor	<a href="#">DWORD</a>	Pen color used to draw trace.
18	nPenWidth	char	Pen width in pixels.
19	bDrawPoints	char	TRUE = Draw disconnected points.
20	bHidden	char	TRUE = Hide the trace.
21	bFloatData	char	TRUE = Floating point pseudo channel.
22	fVertProportion	float	Relative proportion of the client are to use.
26	fDisplayGain	char	Display gain of trace in UserUnits.
30	fDisplayOffset	char	Display offset of trace in UserUnits.
<u>TOTAL 34 bytes</u>			

## # \$ K + The ABF Data Section

Acquired data samples are stored as multiplexed two-byte binary integers. A special four-byte floating point version is used by some analysis programs for storage of analysis results (see [nDataFormat](#) in the file header section). If the data file contains multiple [sweeps](#) or segments of data, these are stored end for end without a gap. That is, no parameters are stored between sweeps of the data; all parameters are stored in the header or in the specialized sections described below. There is only one data section.

---

# The\_ABF\_Data\_Section

\$ The ABF Data Section

K The ABF Data Section

+ ABFFUNC

## # \$ K + The ABF Synch Section

The ABF Synch array is an important array that stores the start time and length of each portion of the data if the data are not part of a continuous [gap-free](#) acquisition. The data section might contain equal length or variable length [sweeps](#) of data. The Synch Array contains a record to indicate the start time and length of every sweep or Event in the data [file](#). The ABF reading routines automatically decode the Synch Array when providing information about the data.

A Synch array is created and used in the following acquisition modes: ABF\_VARLENEVENTS, ABF\_FIXLENEVENTS & ABF\_HIGHSPEEDOSC. The acquisition modes ABF\_GAPFREEFILE and ABF\_WAVEFORMFILE do not always use a Synch array.

### The ABFSynch Structure

(8 bytes).

Offset	Header Entry Name	Type	Description
0	lStart	<a href="#">long</a>	Start time of <a href="#">sweep</a> in fSynchTimeUnit units.
4	lLength	<a href="#">long</a>	Length of the sweep in multiplexed samples.

---

# The\_ABF\_Synch\_Section

\$ The ABF Synch Section

K The ABF Synch Section;ABFSynch

+ ABFFUNC

## # \$ K + The ABF Tag Section

During an acquisition, some programs allow the user to tag points of interest in the input data stream. These tags are saved in the Tag Section. Each tag consists of a time stamp, a text comment, and a tag type identifier. If the tag is a voice tag, the data is held in an ABFVoiceTagInfo struct.

### The ABFTag Structure

(64 bytes).

Offset	Header Entry Name	Type	Description
0	lTagTime	<a href="#">long</a>	Time at which the tag was entered in fSynchTimeUnit units.
4	sComment	56 <a href="#">char</a>	Optional comment to describe the tag.
60	nTagType	<a href="#">short</a>	Type of tag. Valid types are ABF_TIMETAG=0, ABF_COMMENTTAG=1, ABF_EXTERNALTAG=2, ABF_VOICETAG=3
62	nVoiceTagNumber	<a href="#">short</a>	If nTagType=ABF_VOICETAG, this is the number of this voice tag.

### The ABFVoiceTagInfo structure

(32 bytes).

Offset	Header Entry Name	Type	Description
0	lTagNumber	<a href="#">long</a>	The tag number that corresponds to this VoiceTag
4	lFileOffset	long	Offset to this tag within the VoiceTag block
8	lUncompressedSize	long	Size of the voice tag expanded.
12	lCompressedSize	long	Compressed size of the tag.
16	nCompressionType	<a href="#">short</a>	Compression method used.
18	nSampleSize	short	Size of the samples acquired.
20	lSamplesPerSecond	long	Rate at which the sound was acquired.
24	dwCRC	<a href="#">DWORD</a>	CRC used to check data integrity.
28	wChannels	<a href="#">WORD</a>	Number of channels in the tag (usually 1).
30	wUnused	WORD	Unused space.

---

# The\_ABF\_Tag\_Section

\$ The ABF Tag Section

K The ABF Tag Section;ABFTag

+ ABFFUNC

## # \$ K + The ABF Deltas Section

When acquisition parameters are changed during an acquisition, the changes are tracked and entered in the ABF deltas section. Each entry is time stamped in fSynchTimeUnit units, so that the value of the parameter can be determined at any point during the acquisition.

### The ABFDelta Structure

(64 bytes).

Offset	Header Entry Name	Type	Description
0	IDeltaTime	<a href="#">long</a>	Time at which the parameter was changed in fSynchTimeUnit units.
4	IParameterID	<a href="#">long</a>	Identifier for the parameter changed. Legal parameter values are: ABF_DELTA_XXXXXXX
8	INewParamValue fNewParamValue	<a href="#">long</a> <a href="#">float</a>	Depending on the value of IParameterID this entry may be either a <a href="#">float</a> or a <a href="#">long</a> .

---

# The\_ABF\_Deltas\_Section

\$ The ABF Deltas Section

K The ABF Deltas Section;ABFDelta

+ ABFFUNC

## # \$ K + The DAC Data Section

In some experiments, instead of the analog output stimulation waveform being described parametrically in a table it is described on a [sample](#)-by-sample basis from a file. In experiments where files are used to describe the analog output stimulation, a copy of the stimulation file is attached in the [DAC](#) Data section to ensure that the stimulus waveform is always available during analysis. The ABF routines allow reading of the [DAC](#) data with separate functions similar to those of [ABF\\_ReadChannel](#) (i.e. called [ABF\\_ReadDACFileEpi](#)).

---

# The\_DAC\_Data\_Section

\$ The DAC Data Section

K The DAC Data Section

+ ABFFUNC

## # \$ K + ABFHEADR.H

```
//*****
//
// Copyright (c) 1993-2003 Axon Instruments.
// All rights reserved.
// Permission is granted to freely use, modify and copy the code in this file.
//
//*****
// HEADER: ABFHEADR.H.
// PURPOSE: Defines the ABFFileHeader structure, and provides prototypes for
// functions implemented in ABFHEADR.CPP for reading and writing
// ABFFileHeader's.
// REVISIONS:
// 1.1 - Version 1.1 was released in April 1992.
// 1.2 - Added nDataFormat so that data can optionally be stored in floating point format.
// - Added lClockChange to control the multiplexed ADC sample number after which the second
sampling interval commences.
// 1.3 - Change 4-byte sFileType string to long lFileSignature.
// - #define ABF_NATIVESIGNATURE & ABF_REVERSESIGNATURE for byte order detection.
// - Added support for Bells during before or after acquisitions
// - Added parameters to describe hysteresis during event detected acquisitions:
nLevelHysteresis and lTimeHysteresis.
// - Dropped support for BASIC and Pascal.
// - Added the ABF Scope Config section to store scope configuration information
// 1.4 - Remove support for big-endian machines.
// 1.5 - Change ABFSignal parameters from UUTop & UUBottom to
fDisplayGain & fDisplayOffset.
// - Added and changed parameters in the 'File Structure', 'Display Parameters',
// 'DAC Output File', 'Autopeak Measurements' and 'Unused space and end of header' sections
of the ABF file header.
// - Expanded the ABF API and error return codes
// 1.6 - Expanded header to 5120 bytes and added extra parameters to support 2 waveform channels
PRC
// 1.65 - Telegraph support added.
// 1.67 - Train epochs, multiple channel and multiple region stats
// 1.68 - ABFScopeConfig expanded
// 1.69 - Added user entered percentile levels for rise and decay stats
// 1.70 - Added data reduction - AjD
// 1.71 - Added epoch resistance
// 1.72 - Added alternating outputs
// 1.73 - Added post-processing lowpass filter settings. When filtering is done in Clampfit it is
stored in the header.
// 1.74 - Added channel_count_acquired
// 1.75 - Added polarity for each channel
// 1.76 - Added digital trigger out flag
// 1.77 - Added major, minor and bugfix version numbers
// 1.78 - Added separate entries for alternating DAC and digital outputs
// 1.79 - Removed data reduction (now minidigi only)
// 1.80 - Added stats mode for each region: mode is cursor region, epoch etc

#ifndef INC_ABFHEADR_H
#define INC_ABFHEADR_H

#include "AxABFFIO32.h"

#ifdef __cplusplus
extern "C" {
#endif

//
// Constants used in defining the ABF file header
//

#define ABF_ADCCOUNT 16 // number of ADC channels supported.
#define ABF_DACCOUNT 4 // number of DAC channels supported.
#define ABF_WAVEFORMCOUNT 2 // number of DAC channels which support waveforms.
#define ABF_EPOCHCOUNT 10 // number of waveform epochs supported.
```

---

# ABFHEADR\_H

\$ ABFHEADR.H

K ABFHEADR.H

+ ABFFUNC



```

#define ABF_BELLCOUNT          2      // Number of auditory signals supported.
#define ABF_ADCUNITLEN         8      // length of ADC units strings
#define ABF_ADCNAMELEN         10     // length of ADC channel name strings
#define ABF_DACUNITLEN         8      // length of DAC units strings
#define ABF_DACNAMELEN         10     // length of DAC channel name strings
#define ABF_VARPARAMLISTLEN    80     // length of conditioning string
#define ABF_USERLISTLEN        256    // length of the user list (V1.6)
#define ABF_USERLISTCOUNT     4      // number of independent user lists (V1.6)
#define ABF_OLDFILECOMMENTLEN  56     // length of file comment string (pre V1.6)
#define ABF_FILECOMMENTLEN     128    // length of file comment string (V1.6)

#define ABF_CREATORINFOLEN      16     // length of file creator info string
#define ABF_OLDDACFILENAMELEN   12     // old length of the DACFile name string
#define ABF_OLDDACFILEPATHLEN  60     // old length of the DACFile path string
#define ABF_DACFILEPATHLEN     84     // length of full path for DACFile
#define ABF_PATHLEN             256    // length of full path, used for DACFile and Protocol name.
#define ABF_ARITHMETICOPLEN     2      // length of the Arithmetic operator field
#define ABF_ARITHMETICUNITSLEN  8      // length of arithmetic units string
#define ABF_TAGCOMMENTLEN       56     // length of tag comment string
#define ABF_LONGDESCRIPTIONLEN  56     // length of long description entry
#define ABF_NOTENAMELEN         10     // length of the name component of a note
#define ABF_NOTEVALUELEN        8      // length of the value component of a note
#define ABF_NOTEUNITSLEN        8      // length of the units component of a note
#define ABF_BLOCKSIZE           512    // Size of block alignment in ABF files.
#define ABF_MACRONAMELEN        64     // Size of a Clampfit macro name.

#define ABF_CURRENTVERSION      1.80F  // Current file format version number
#define ABF_PREVIOUSVERSION     1.5F   // Previous file format version number (for old
header size)
#define ABF_V16                 1.6F   // Version number when the header size changed.
#define ABF_HEADERSIZE          6144   // Size of a Version 1.6 or later header
#define ABF_OLDHEADERSIZE       2048   // Size of a Version 1.5 or earlier header
#define ABF_NATIVESIGNATURE     0x20464241 // PC="ABF ", MAC=" FBA"
#define ABF_REVERSESIGNATURE    0x41424620 // PC=" FBA", MAC="ABF "

#define PCLAMP6_MAXSWEEPLENGTH  16384  // Maximum multiplexed sweep length supported by
pCLAMP6 apps.
#define PCLAMP7_MAXSWEEPLEN_PERCHAN 1032258 // Maximum per channel sweep length supported by
pCLAMP7 apps.

#define ABF_MAX_TRIAL_SAMPLES  0x7FFFFFFF // Maximum length of acquisition supported (samples)
// INT_MAX is used instead of UINT_MAX because of the
signed
// values in the ABF header.

#define ABF_MAX_SWEEPS_PER_AVERAGE 65500 // The maximum number of sweeps that can be combined
into a
// cumulative average
(nAverageAlgorithm=ABF_INFINITEAVERAGE).

#define ABF_STATS_REGIONS      8      // The number of independent statistics regions.
#define ABF_BASELINE_REGIONS   1      // The number of independent baseline regions.

#ifdef _MAC
#define ABF_OLDPCLAMP          ABF_REVERSESIGNATURE
#else
#define ABF_OLDPCLAMP          ABF_NATIVESIGNATURE
#endif

//
// Constant definitions for nFileType
//
#define ABF_ABFFILE             1
#define ABF_FETCHEX             2
#define ABF_CLAMPEX             3

//
// Constant definitions for nDataFormat
//
#define ABF_INTEGERDATA         0
#define ABF_FLOATDATA           1

//
// Constant definitions for nOperationMode
//
#define ABF_VARLENEVENTS        1
#define ABF_FIXLENEVENTS        2      // (ABF_FIXLENEVENTS == ABF_LOSSFREEOSC)
#define ABF_LOSSFREEOSC         2
#define ABF_GAPFREEFILE          3

```

```

#define ABF_HIGHSPEEDOSC      4
#define ABF_WAVEFORMFILE     5

//
// Constant definitions for nParamToVary
//
#define ABF_CONDITNUMPULSES      0
#define ABF_CONDITBASELINEDURATION 1
#define ABF_CONDITBASELINELEVEL 2
#define ABF_CONDITSTEPDURATION  3
#define ABF_CONDITSTEPLEVEL     4
#define ABF_CONDITPOSTTRAINDURATION 5
#define ABF_CONDITPOSTTRAINLEVEL 6
#define ABF_EPISODESTARTTOSTART 7
#define ABF_INACTIVEHOLDING      8
#define ABF_DIGITALHOLDING       9
#define ABF_PNNUMPULSES         10
#define ABF_PARALLELVALUE        11
#define ABF_EPOCHINITLEVEL      (ABF_PARALLELVALUE + ABF_EPOCHCOUNT)
#define ABF_EPOCHINITDURATION   (ABF_EPOCHINITLEVEL + ABF_EPOCHCOUNT)
#define ABF_EPOCHTRAINPERIOD    (ABF_EPOCHINITDURATION + ABF_EPOCHCOUNT)
#define ABF_EPOCHTRAINPULSEWIDTH (ABF_EPOCHTRAINPERIOD + ABF_EPOCHCOUNT)
// Next value is (ABF_EPOCHINITDURATION + ABF_EPOCHCOUNT)

//
// Constants for nAveragingMode
//
#define ABF_NOAVERAGING          0
#define ABF_SAVEAVERAGEONLY    1
#define ABF_AVERAGESAVEALL      2

//
// Constants for nAverageAlgorithm
//
#define ABF_INFINITEAVERAGE     0
#define ABF_SLIDINGAVERAGE      1

//
// Constants for nEpochType
//
#define ABF_EPOCHDISABLED        0    // disabled epoch
#define ABF_EPOCHSTEPPED         1    // stepped waveform
#define ABF_EPOCHRAMPED          2    // ramp waveform
#define ABF_EPOCH_TYPE_RECTANGLE 3    // rectangular pulse train
#define ABF_EPOCH_TYPE_TRIANGLE 4    // triangular waveform
#define ABF_EPOCH_TYPE_COSINE    5    // cosinusoidal waveform
#define ABF_EPOCH_TYPE_RESISTANCE 6    // resistance waveform
#define ABF_EPOCH_TYPE_BIPHASIC  7    // biphasic pulse train

//
// Constants for epoch resistance
//
#define ABF_MIN_EPOCH_RESISTANCE_DURATION 8

//
// Constants for nWaveformSource
//
#define ABF_WAVEFORMDISABLED      0    // disabled waveform
#define ABF_EPOCHTABLEWAVEFORM   1
#define ABF_DACFILEWAVEFORM      2

//
// Constants for nInterEpisodeLevel & nDigitalInterEpisode
//
#define ABF_INTEREPI_USEHOLDING   0
#define ABF_INTEREPI_USELASTEPOCH 1

//
// Constants for nExperimentType
//
#define ABF_VOLTAGECLAMP          0
#define ABF_CURRENTCLAMP         1
#define ABF_SIMPLEACQUISITION    2

//
// Constants for nAutosampleEnable
//
#define ABF_AUTOSAMPLEDISABLED    0
#define ABF_AUTOSAMPLEAUTOMATIC  1

```

```

#define ABF_AUTOSAMPLEMANUAL      2

//
// Constants for nAutosampleInstrument
//
#define ABF_INST_UNKNOWN          0 // Unknown instrument (manual or user defined telegraph table).
#define ABF_INST_AXOPATCH1        1 // Axopatch-1 with CV-4-1/100
#define ABF_INST_AXOPATCH1_1      2 // Axopatch-1 with CV-4-0.1/100
#define ABF_INST_AXOPATCH1B       3 // Axopatch-1B(inv.) CV-4-1/100
#define ABF_INST_AXOPATCH1B_1     4 // Axopatch-1B(inv) CV-4-0.1/100
#define ABF_INST_AXOPATCH201      5 // Axopatch 200 with CV 201
#define ABF_INST_AXOPATCH202      6 // Axopatch 200 with CV 202
#define ABF_INST_GENECLAMP        7 // GeneClamp
#define ABF_INST_DAGAN3900        8 // Dagan 3900
#define ABF_INST_DAGAN3900A       9 // Dagan 3900A
#define ABF_INST_DAGANCA1_1       10 // Dagan CA-1 Im=0.1
#define ABF_INST_DAGANCA1        11 // Dagan CA-1 Im=1.0
#define ABF_INST_DAGANCA10       12 // Dagan CA-1 Im=10
#define ABF_INST_WARNER_OC725     13 // Warner OC-725
#define ABF_INST_WARNER_OC725C    14 // Warner OC-725
#define ABF_INST_AXOPATCH200B     15 // Axopatch 200B
#define ABF_INST_DAGANPCONE0_1     16 // Dagan PC-ONE Im=0.1
#define ABF_INST_DAGANPCONE1      17 // Dagan PC-ONE Im=1.0
#define ABF_INST_DAGANPCONE10     18 // Dagan PC-ONE Im=10
#define ABF_INST_DAGANPCONE100    19 // Dagan PC-ONE Im=100
#define ABF_INST_WARNER_BC525C    20 // Warner BC-525C
#define ABF_INST_WARNER_PC505     21 // Warner PC-505
#define ABF_INST_WARNER_PC501     22 // Warner PC-501
#define ABF_INST_DAGANCA1_05      23 // Dagan CA-1 Im=0.05
#define ABF_INST_MULTICLAMP700    24 // MultiClamp 700
#define ABF_INST_TURBO_TEC        25 // Turbo Tec
#define ABF_INST_OPUSXPRESS6000   26 // OpusXpress 6000A

//
// Constants for nManualInfoStrategy
//
#define ABF_ENV_DONOTWRITE        0
#define ABF_ENV_WRITEEACHTRIAL    1
#define ABF_ENV_PROMPTTEACHTRIAL  2

//
// Constants for nTriggerSource
//
#define ABF_TRIGGERLINEINPUT      -5 // Start on line trigger (DD1320 only)
#define ABF_TRIGGERTAGINPUT       -4
#define ABF_TRIGGERFIRSTCHANNEL   -3
#define ABF_TRIGGEREXTERNAL       -2
#define ABF_TRIGGERSPACEBAR       -1
// >=0 = ADC channel to trigger off.

//
// Constants for nTrialTriggerSource
//
#define ABF_TRIALTRIGGER_SWSTARTONLY -6 // Start on software message, end when protocol ends.
#define ABF_TRIALTRIGGER_SWSTARTSTOP -5 // Start and end on software messages.
#define ABF_TRIALTRIGGER_LINEINPUT  -4 // Start on line trigger (DD1320 only)
#define ABF_TRIALTRIGGER_SPACEBAR    -3 // Start on spacebar press.
#define ABF_TRIALTRIGGER_EXTERNAL    -2 // Start on external trigger high
#define ABF_TRIALTRIGGER_NONE        -1 // Start immediately (default).
// >=0 = ADC channel to trigger off. // Not implemented as yet...

//
// Constants for nTriggerPolarity.
//
#define ABF_TRIGGER_RISINGEDGE      0
#define ABF_TRIGGER_FALLINGEDGE     1

//
// Constants for nTriggerAction
//
#define ABF_TRIGGER_STARTEPISODE    0
#define ABF_TRIGGER_STARTRUN        1
#define ABF_TRIGGER_STARTTRIAL      2 // N.B. Discontinued in favor of nTrialTriggerSource

//
// Constants for nDrawingStrategy
//
#define ABF_DRAW_NONE                0
#define ABF_DRAW_REALTIME            1

```

```

#define ABF_DRAW_FULLSCREEN      2
#define ABF_DRAW_ENDOFRUN      3

//
// Constants for nTiledDisplay
//
#define ABF_DISPLAY_SUPERIMPOSED 0
#define ABF_DISPLAY_TILED      1

//
// Constants for nDataDisplayMode
//
#define ABF_DRAW_POINTS        0
#define ABF_DRAW_LINES        1

//
// Constants for nArithmeticExpression
//
#define ABF_SIMPLE_EXPRESSION    0
#define ABF_RATIO_EXPRESSION    1

//
// Constants for nLowpassFilterType & nHighpassFilterType
//
#define ABF_FILTER_NONE          0
#define ABF_FILTER_EXTERNAL      1
#define ABF_FILTER_SIMPLE_RC     2
#define ABF_FILTER_BESSEL       3
#define ABF_FILTER_BUTTERWORTH  4

//
// Constants for nPNPosition
//
#define ABF_PN_BEFORE_EPISODE    0
#define ABF_PN_AFTER_EPISODE    1

//
// Constants for nPNPolarity
//
#define ABF_PN_OPPOSITE_POLARITY -1
#define ABF_PN_SAME_POLARITY     1

//
// Constants for nAutopeakPolarity
//
#define ABF_PEAK_NEGATIVE        -1
#define ABF_PEAK_ABSOLUTE        0
#define ABF_PEAK_POSITIVE        1

//
// Constants for nAutopeakSearchMode
//
#define ABF_PEAK_SEARCH_SPECIFIED -2
#define ABF_PEAK_SEARCH_ALL      -1
// nAutopeakSearchMode 0..9 = epoch in waveform 0's epoch table
// nAutopeakSearchMode 10..19 = epoch in waveform 1's epoch table

//
// Constants for nAutopeakBaseline
//
#define ABF_PEAK_BASELINE_SPECIFIED -3
#define ABF_PEAK_BASELINE_NONE      -2
#define ABF_PEAK_BASELINE_FIRSTHOLDING -1
#define ABF_PEAK_BASELINE_LASTHOLDING -4

//
// Constants for lAutopeakMeasurements
//
#define ABF_PEAK_MEASURE_PEAK          0x00000001
#define ABF_PEAK_MEASURE_PEAKTIME      0x00000002
#define ABF_PEAK_MEASURE_ANTIPEAK      0x00000004
#define ABF_PEAK_MEASURE_ANTIPEAKTIME  0x00000008
#define ABF_PEAK_MEASURE_MEAN          0x00000010
#define ABF_PEAK_MEASURE_STDDEV        0x00000020
#define ABF_PEAK_MEASURE_INTEGRAL      0x00000040
#define ABF_PEAK_MEASURE_MAXRISESLOPE  0x00000080
#define ABF_PEAK_MEASURE_MAXRISESLOPETIME 0x00000100
#define ABF_PEAK_MEASURE_MAXDECAYSLOPE 0x00000200
#define ABF_PEAK_MEASURE_MAXDECAYSLOPETIME 0x00000400

```

```

#define ABF_PEAK_MEASURE_RISETIME          0x00000800
#define ABF_PEAK_MEASURE_DECAYTIME         0x00001000
#define ABF_PEAK_MEASURE_HALFWIDTH        0x00002000
#define ABF_PEAK_MEASURE_BASELINE          0x00004000
#define ABF_PEAK_MEASURE_RISESLOPE         0x00008000
#define ABF_PEAK_MEASURE_DECAYSLOPE        0x00010000
#define ABF_PEAK_MEASURE_REGIONSLOPE       0x00020000
#define ABF_PEAK_MEASURE_ALL                0x0002FFFF // All of the above OR'd together.

//
// Constants for nStatsActiveChannels
//
#define ABF_PEAK_SEARCH_CHANNEL0           0x0001
#define ABF_PEAK_SEARCH_CHANNEL1           0x0002
#define ABF_PEAK_SEARCH_CHANNEL2           0x0004
#define ABF_PEAK_SEARCH_CHANNEL3           0x0008
#define ABF_PEAK_SEARCH_CHANNEL4           0x0010
#define ABF_PEAK_SEARCH_CHANNEL5           0x0020
#define ABF_PEAK_SEARCH_CHANNEL6           0x0040
#define ABF_PEAK_SEARCH_CHANNEL7           0x0080
#define ABF_PEAK_SEARCH_CHANNEL8           0x0100
#define ABF_PEAK_SEARCH_CHANNEL9           0x0200
#define ABF_PEAK_SEARCH_CHANNEL10          0x0400
#define ABF_PEAK_SEARCH_CHANNEL11          0x0800
#define ABF_PEAK_SEARCH_CHANNEL12          0x1000
#define ABF_PEAK_SEARCH_CHANNEL13          0x2000
#define ABF_PEAK_SEARCH_CHANNEL14          0x4000
#define ABF_PEAK_SEARCH_CHANNEL15          0x8000
#define ABF_PEAK_SEARCH_CHANNELSALL        0xFFFF // All of the above OR'd together.

// Bit flag settings for nStatsSearchRegionFlags
//
#define ABF_PEAK_SEARCH_REGION0            0x01
#define ABF_PEAK_SEARCH_REGION1            0x02
#define ABF_PEAK_SEARCH_REGION2            0x04
#define ABF_PEAK_SEARCH_REGION3            0x08
#define ABF_PEAK_SEARCH_REGION4            0x10
#define ABF_PEAK_SEARCH_REGION5            0x20
#define ABF_PEAK_SEARCH_REGION6            0x40
#define ABF_PEAK_SEARCH_REGION7            0x80
#define ABF_PEAK_SEARCH_REGIONALL          0xFF // All of the above OR'd together.

//
// Constants for lStatisticsMeasurements
//
#define ABF_STATISTICS_ABOVETHRESHOLD       0x00000001
#define ABF_STATISTICS_EVENTFREQUENCY       0x00000002
#define ABF_STATISTICS_MEANOPENTIME         0x00000004
#define ABF_STATISTICS_MEANCLOSEDTIME       0x00000008
#define ABF_STATISTICS_ALL                  0x0000000F // All the above OR'd together.

//
// Constants for nStatisticsSaveStrategy
//
#define ABF_STATISTICS_NOAUTOSAVE            0
#define ABF_STATISTICS_AUTOSAVE              1
#define ABF_STATISTICS_AUTOSAVE_AUTOCLEAR    2

//
// Constants for nStatisticsDisplayStrategy
//
#define ABF_STATISTICS_DISPLAY               0
#define ABF_STATISTICS_NODISPLAY            1

//
// Constants for nStatisticsClearStrategy
// determines whether to clear statistics after saving.
//
#define ABF_STATISTICS_NOCLEAR               0
#define ABF_STATISTICS_CLEAR                 1

//
// Constants for nDACFileEpisodeNum
//
#define ABF_DACFILE_SKIPFIRSTSWEEP -1
#define ABF_DACFILE_USEALLSWEEPS 0
// >0 = The specific sweep number.

//

```

```

// Constants for nUndoPromptStrategy
//
#define ABF_UNDOPROMPT_ONABORT 0
#define ABF_UNDOPROMPT_ALWAYS 1

//
// Constants for nAutoAnalyseEnable
//
#define ABF_AUTOANALYSE_DISABLED 0
#define ABF_AUTOANALYSE_DEFAULT 1
#define ABF_AUTOANALYSE_RUNMACRO 2

//
// Constants for post nPostprocessLowpassFilterType
//
#define ABF_POSTPROCESS_FILTER_NONE 0
#define ABF_POSTPROCESS_FILTER_COMBINATION 1
#define ABF_POSTPROCESS_FILTER_BESSEL 2
#define ABF_POSTPROCESS_FILTER_BOXCAR 3
#define ABF_POSTPROCESS_FILTER_BUTTERWORTH 4
#define ABF_POSTPROCESS_FILTER_CHEBYSHEV 5
#define ABF_POSTPROCESS_FILTER_GAUSSIAN 6
#define ABF_POSTPROCESS_FILTER_RC 7
#define ABF_POSTPROCESS_FILTER_RC8 8
#define ABF_POSTPROCESS_FILTER_ADAPTIVE 9

//
// Miscellaneous constants
//
#define ABF_FILTERDISABLED 100000.0F // Large frequency to disable lowpass filters
#define ABF_UNUSED_CHANNEL -1 // Unused ADC and DAC channels.

//
// The output sampling sequence identifier for a separate digital out channel.
//
#define ABF_DIGITAL_OUT_CHANNEL -1
#define ABF_PADDING_OUT_CHANNEL -2

//
// maximum values for various parameters (used by ABFH_CheckUserList).
//
#define ABF_CTPULSECOUNT_MAX 10000
#define ABF_CTBASELINEDURATION_MAX 100000.0F
#define ABF_CTSTEPDURATION_MAX 100000.0F
#define ABF_CTPOSTTRAINDURATION_MAX 100000.0F
#define ABF_SWEEPSTARTTOSTARTTIME_MAX 100000.0F
#define ABF_PNPULSECOUNT_MAX 8
#define ABF_DIGITALVALUE_MAX 0xFF
#define ABF_EPOCHDIGITALVALUE_MAX 0x0F

//
// LTP Types - Reflects whether the header is used for LTP as baseline or induction.
//
#define ABF_LTP_TYPE_NONE 0
#define ABF_LTP_TYPE_BASELINE 1
#define ABF_LTP_TYPE_INDUCTION 2

//
// LTP Usage of DAC - Reflects whether the analog output will be used presynaptically or
// postsynaptically.
//
#define ABF_LTP_DAC_USAGE_NONE 0
#define ABF_LTP_DAC_USAGE_PRESYNAPTIC 1
#define ABF_LTP_DAC_USAGE_POSTSYNAPTIC 2

//
// Header Version Numbers
//
#define ABF_V166 1.66F
#define ABF_V167 1.67F
#define ABF_V168 1.68F
#define ABF_V169 1.69F
#define ABF_V170 1.70F
#define ABF_V171 1.71F
#define ABF_V172 1.72F
#define ABF_V173 1.73F
#define ABF_V174 1.74F
#define ABF_V175 1.75F
#define ABF_V176 1.76F

```

```

#define ABF_V177 1.77F
#define ABF_V178 1.78F
#define ABF_V179 1.79F
#define ABF_V180 1.80F

//
// pack structure on byte boundaries
//
#ifndef RC_INVOKED
#pragma pack(push, 1)
#endif

//
// Definition of the ABF header structure.
//

struct ABFFileHeader // The total header length = 6144 bytes.
{
public:
    // GROUP #1 - File ID and size information. (40 bytes)
    long    lFileSignature;
    float    fFileVersionNumber;
    short    nOperationMode;
    long    lActualAcqLength;
    short    nNumPointsIgnored;
    long    lActualEpisodes;
    long    lFileStartDate; // YYYYMMDD
    long    lFileStartTime;
    long    lStopwatchTime;
    float    fHeaderVersionNumber;
    short    nFileType;
    short    nMSBinFormat;

    // GROUP #2 - File Structure (78 bytes)
    long    lDataSectionPtr;
    long    lTagSectionPtr;
    long    lNumTagEntries;
    long    lScopeConfigPtr;
    long    lNumScopes;
    long    _lDACFilePtr;
    long    _lDACFileNumEpisodes;
    char    sUnused001[4];
    long    lDeltaArrayPtr;
    long    lNumDeltas;
    long    lVoiceTagPtr;
    long    lVoiceTagEntries;
    long    lUnused002;
    long    lSynchArrayPtr;
    long    lSynchArraySize;
    short    nDataFormat;
    short    nSimultaneousScan;
    long    lStatisticsConfigPtr;
    long    lAnnotationSectionPtr;
    long    lNumAnnotations;
    char    sUnused003[2];

    // GROUP #3 - Trial hierarchy information (82 bytes)
    /**
    The number of input channels we acquired.
    Do not access directly - use CABFHeader::get_channel_count_acquired
    */
    short    channel_count_acquired;

    /**
    The number of input channels we recorded.
    Do not access directly - use CABFHeader::get_channel_count_recorded
    */
    short    nADCNumChannels;
    float    fADCSampleInterval;
    /*{
    The documentation says these two sample intervals are the interval between multiplexed
    samples, but not all digitisers work like that.
    Instead, these are the per-channel sample rate divided by the number of channels.
    If the user chose 100uS and has two channels, this value will be 50uS.
    }*/
    float    fADCSecondSampleInterval;
    /*{
    // The two sample intervals must be an integer multiple (or submultiple) of each other.
    if (fADCSampleInterval > fADCSecondSampleInterval)

```

```

        ASSERT(fmod(fADCSampleInterval, fADCSecondSampleInterval) == 0.0);
    if (fADCSecondSampleInterval, fADCSampleInterval)
        ASSERT(fmod(fADCSecondSampleInterval, fADCSampleInterval) == 0.0);
    }*/
float    fSynchTimeUnit;
float    fSecondsPerRun;

/**
 * The total number of samples per episode, for the recorded channels only.
 * This does not include channels which are acquired but not recorded.
 *
 * This is the number of samples per episode per channel, times the number of recorded channels.
 *
 * If you want the samples per episode for one channel, you must divide this by
    get_channel_count_recorded().
 */
long      lNumSamplesPerEpisode;
long      lPreTriggerSamples;
long      lEpisodesPerRun;
long      lRunsPerTrial;
long      lNumberOfTrials;
short     nAveragingMode;
short     nUndoRunCount;
short     nFirstEpisodeInRun;
float     fTriggerThreshold;
short     nTriggerSource;
short     nTriggerAction;
short     nTriggerPolarity;
float     fScopeOutputInterval;
float     fEpisodeStartToStart;
float     fRunStartToStart;
float     fTrialStartToStart;
long      lAverageCount;
long      lClockChange;
short     nAutoTriggerStrategy;

// GROUP #4 - Display Parameters (44 bytes)
short     nDrawingStrategy;
short     nTiledDisplay;
short     nEraseStrategy;           // N.B. Discontinued. Use scope config entry instead.
short     nDataDisplayMode;
long      lDisplayAverageUpdate;
short     nChannelStatsStrategy;
long      lCalculationPeriod;       // N.B. Discontinued. Use fStatisticsPeriod.
long      lSamplesPerTrace;
long      lStartDisplayNum;
long      lFinishDisplayNum;
short     nMultiColor;
short     nShowPNRawData;
float     fStatisticsPeriod;
long      lStatisticsMeasurements;
short     nStatisticsSaveStrategy;

// GROUP #5 - Hardware information (16 bytes)
float     fADCRange;
float     fDACRange;
long      lADCResolution;
long      lDACResolution;

// GROUP #6 Environmental Information (118 bytes)
short     nExperimentType;
short     _nAutosampleEnable;
short     _nAutosampleADCNum;
short     _nAutosampleInstrument;
float     _fAutosampleAdditGain;
float     _fAutosampleFilter;
float     _fAutosampleMembraneCap;
short     nManualInfoStrategy;
float     fCellID1;
float     fCellID2;
float     fCellID3;
char      sCreatorInfo[ABF_CREATORINFOLEN];
char      _sFileComment[ABF_OLDFILECOMMENTLEN];
short     nFileStartMillisecs;     // Milliseconds portion of lFileStartTime
short     nCommentsEnable;
char      sUnused003a[8];

// GROUP #7 - Multi-channel information (1044 bytes)
short     nADCPTolChannelMap[ABF_ADCCOUNT];

```



```

short    nADCSamplingSeq[ABF_ADCCOUNT];
char     sADCChannelName[ABF_ADCCOUNT][ABF_ADCNAMELEN];
char     sADCUnits[ABF_ADCCOUNT][ABF_ADUNITLEN];
float    fADCProgrammableGain[ABF_ADCCOUNT];
float    fADCDisplayAmplification[ABF_ADCCOUNT];
float    fADCDisplayOffset[ABF_ADCCOUNT];
float    fInstrumentScaleFactor[ABF_ADCCOUNT];
float    fInstrumentOffset[ABF_ADCCOUNT];
float    fSignalGain[ABF_ADCCOUNT];
float    fSignalOffset[ABF_ADCCOUNT];
float    fSignalLowpassFilter[ABF_ADCCOUNT];
float    fSignalHighpassFilter[ABF_ADCCOUNT];
char     sDACChannelName[ABF_DACCOUNT][ABF_DACNAMELEN];
char     sDACChannelUnits[ABF_DACCOUNT][ABF_DACUNITLEN];
float    fDACScaleFactor[ABF_DACCOUNT];
float    fDACHoldingLevel[ABF_DACCOUNT];
short    nSignalType;
char     sUnused004[10];

// GROUP #8 - Synchronous timer outputs (14 bytes)
short    nOUTEnable;
short    nSampleNumberOUT1;
short    nSampleNumberOUT2;
short    nFirstEpisodeOUT;
short    nLastEpisodeOUT;
short    nPulseSamplesOUT1;
short    nPulseSamplesOUT2;

// GROUP #9 - Epoch Waveform and Pulses (184 bytes)
short    nDigitalEnable;
short    _nWaveformSource;
short    nActiveDACChannel;
short    _nInterEpisodeLevel;
short    _nEpochType[ABF_EPOCHCOUNT];
float    _fEpochInitLevel[ABF_EPOCHCOUNT];
float    _fEpochLevelInc[ABF_EPOCHCOUNT];
short    _nEpochInitDuration[ABF_EPOCHCOUNT];
short    _nEpochDurationInc[ABF_EPOCHCOUNT];
short    nDigitalHolding;
short    nDigitalInterEpisode;
short    nDigitalValue[ABF_EPOCHCOUNT];
char     sUnavailable1608[4];    // was float fWaveformOffset;
short    nDigitalDACChannel;
char     sUnused005[6];

// GROUP #10 - DAC Output File (98 bytes)
float    _fDACFileScale;
float    _fDACFileOffset;
char     sUnused006[2];
short    _nDACFileEpisodeNum;
short    _nDACFileADCNum;
char     _sDACFilePath[ABF_DACFILEPATHLEN];

// GROUP #11 - Presweep (conditioning) pulse train (44 bytes)
short    _nConditEnable;
short    _nConditChannel;
long     _lConditNumPulses;
float    _fBaselineDuration;
float    _fBaselineLevel;
float    _fStepDuration;
float    _fStepLevel;
float    _fPostTrainPeriod;
float    _fPostTrainLevel;
char     sUnused007[12];

// GROUP #12 - Variable parameter user list ( 82 bytes)
short    _nParamToVary;
char     _sParamValueList[ABF_VARPARAMLISTLEN];

// GROUP #13 - Autopeak measurement (36 bytes)
short    _nAutopeakEnable;
short    _nAutopeakPolarity;
short    _nAutopeakADCNum;
short    _nAutopeakSearchMode;
long     _lAutopeakStart;
long     _lAutopeakEnd;
short    _nAutopeakSmoothing;
short    _nAutopeakBaseline;
short    _nAutopeakAverage;

```

```

char    sUnavailable1866[2];    // Was nAutopeakSaveStrategy, use nStatisticsSaveStrategy
long    _lAutopeakBaselineStart;
long    _lAutopeakBaselineEnd;
long    _lAutopeakMeasurements;

// GROUP #14 - Channel Arithmetic (52 bytes)
short    nArithmeticEnable;
float    fArithmeticUpperLimit;
float    fArithmeticLowerLimit;
short    nArithmeticADCNumA;
short    nArithmeticADCNumB;
float    fArithmeticK1;
float    fArithmeticK2;
float    fArithmeticK3;
float    fArithmeticK4;
char    sArithmeticOperator[ABF_ARITHMETICOPLEN];
char    sArithmeticUnits[ABF_ARITHMETICUNITSLEN];
float    fArithmeticK5;
float    fArithmeticK6;
short    nArithmeticExpression;
char    sUnused008[2];

// GROUP #15 - On-line subtraction (34 bytes)
short    _nPNEnable;
short    nPNPosition;
short    _nPNPolarity;
short    nPNNumPulses;
short    _nPNADCNum;
float    _fPNHoldingLevel;
float    fPNSettlingTime;
float    fPNInterpulse;
char    sUnused009[12];

// GROUP #16 - Miscellaneous variables (82 bytes)
short    _nListEnable;

short    nBellEnable[ABF_BELLCOUNT];
short    nBellLocation[ABF_BELLCOUNT];
short    nBellRepetitions[ABF_BELLCOUNT];

short    nLevelHysteresis;
long    lTimeHysteresis;
short    nAllowExternalTags;

char    nLowpassFilterType[ABF_ADCCOUNT];
char    nHighpassFilterType[ABF_ADCCOUNT];
short    nAverageAlgorithm;
float    fAverageWeighting;
short    nUndoPromptStrategy;
short    nTrialTriggerSource;
short    nStatisticsDisplayStrategy;
short    nExternalTagType;
long    lHeaderSize;
double   dFileDuration;
short    nStatisticsClearStrategy;
// Size of v1.5 header = 2048

// Extra parameters in v1.6
// EXTENDED GROUP #2 - File Structure (26 bytes)
long    lDACFilePtr[ABF_WAVEFORMCOUNT];
long    lDACFileNumEpisodes[ABF_WAVEFORMCOUNT];
char    sUnused010[10];

// EXTENDED GROUP #7 - Multi-channel information (62 bytes)
float    fDACCcalibrationFactor[ABF_DACCOUNT];
float    fDACCcalibrationOffset[ABF_DACCOUNT];
char    sUnused011[30];

// GROUP #17 - Trains parameters (160 bytes)
long    lEpochPulsePeriod[ABF_WAVEFORMCOUNT][ABF_EPOCHCOUNT];
long    lEpochPulseWidth [ABF_WAVEFORMCOUNT][ABF_EPOCHCOUNT];

// EXTENDED GROUP #9 - Epoch Waveform and Pulses ( 412 bytes)
short    nWaveformEnable[ABF_WAVEFORMCOUNT];
short    nWaveformSource[ABF_WAVEFORMCOUNT];
short    nInterEpisodeLevel[ABF_WAVEFORMCOUNT];
short    nEpochType[ABF_WAVEFORMCOUNT][ABF_EPOCHCOUNT];
float    fEpochInitLevel[ABF_WAVEFORMCOUNT][ABF_EPOCHCOUNT];
float    fEpochLevelInc[ABF_WAVEFORMCOUNT][ABF_EPOCHCOUNT];

```

```

long    lEpochInitDuration[ABF_WAVEFORMCOUNT][ABF_EPOCHCOUNT];
long    lEpochDurationInc[ABF_WAVEFORMCOUNT][ABF_EPOCHCOUNT];
short   nDigitalTrainValue[ABF_EPOCHCOUNT];           // 2 * 10 = 20 bytes
short   nDigitalTrainActiveLogic;                       // 2 bytes
char     sUnused012[18];

// EXTENDED GROUP #10 - DAC Output File (552 bytes)
float    fDACFileScale[ABF_WAVEFORMCOUNT];
float    fDACFileOffset[ABF_WAVEFORMCOUNT];
long     lDACFileEpisodeNum[ABF_WAVEFORMCOUNT];
short    nDACFileADCNum[ABF_WAVEFORMCOUNT];
char     sDACFilePath[ABF_WAVEFORMCOUNT][ABF_PATHLEN];
char     sUnused013[12];

// EXTENDED GROUP #11 - Presweep (conditioning) pulse train (100 bytes)
short    nConditEnable[ABF_WAVEFORMCOUNT];
long     lConditNumPulses[ABF_WAVEFORMCOUNT];
float    fBaselineDuration[ABF_WAVEFORMCOUNT];
float    fBaselineLevel[ABF_WAVEFORMCOUNT];
float    fStepDuration[ABF_WAVEFORMCOUNT];
float    fStepLevel[ABF_WAVEFORMCOUNT];
float    fPostTrainPeriod[ABF_WAVEFORMCOUNT];
float    fPostTrainLevel[ABF_WAVEFORMCOUNT];
char     sUnused014[40];

// EXTENDED GROUP #12 - Variable parameter user list (1096 bytes)
short    nULEnable[ABF_USERLISTCOUNT];
short    nULParamToVary[ABF_USERLISTCOUNT];
char     sULParamValueList[ABF_USERLISTCOUNT][ABF_USERLISTLEN];
short    nULRepeat[ABF_USERLISTCOUNT];
char     sUnused015[48];

// EXTENDED GROUP #15 - On-line subtraction (56 bytes)
short    nPNEnable[ABF_WAVEFORMCOUNT];
short    nPNPolarity[ABF_WAVEFORMCOUNT];
short    nPNADCNum[ABF_WAVEFORMCOUNT];
float    fPNHoldingLevel[ABF_WAVEFORMCOUNT];
char     sUnused016[36];

// EXTENDED GROUP #6 Environmental Information (898 bytes)
short    nTelegraphEnable[ABF_ADCCOUNT];
short    nTelegraphInstrument[ABF_ADCCOUNT];
float    fTelegraphAdditGain[ABF_ADCCOUNT];
float    fTelegraphFilter[ABF_ADCCOUNT];
float    fTelegraphMembraneCap[ABF_ADCCOUNT];
short    nTelegraphMode[ABF_ADCCOUNT];
short    nTelegraphDACScaleFactorEnable[ABF_DACCOUNT];
char     sUnused016a[24];

short    nAutoAnalyseEnable;
char     sAutoAnalysisMacroName[ABF_MACRONAMELEN];
char     sProtocolPath[ABF_PATHLEN];

char     sFileComment[ABF_FILECOMMENTLEN];
char     sUnused017[128];

// EXTENDED GROUP #13 - Statistics measurements (388 bytes)
short    nStatsEnable;
unsigned short nStatsActiveChannels;           // Active stats channel bit flag
unsigned short nStatsSearchRegionFlags;        // Active stats region bit flag
short    nStatsSelectedRegion;
short    _nStatsSearchMode;
short    nStatsSmoothing;
short    nStatsSmoothingEnable;
short    nStatsBaseline;
long     lStatsBaselineStart;
long     lStatsBaselineEnd;
long     lStatsMeasurements[ABF_STATS_REGIONS]; // Measurement bit flag for each region
long     lStatsStart[ABF_STATS_REGIONS];
long     lStatsEnd[ABF_STATS_REGIONS];
short    nRiseBottomPercentile[ABF_STATS_REGIONS];
short    nRiseTopPercentile[ABF_STATS_REGIONS];
short    nDecayBottomPercentile[ABF_STATS_REGIONS];
short    nDecayTopPercentile[ABF_STATS_REGIONS];
short    nStatsChannelPolarity[ABF_ADCCOUNT];
short    nStatsSearchMode[ABF_STATS_REGIONS];    // Stats mode per region: mode is cursor region,
epoch etc
char     sUnused018[156];

```

```

// GROUP #18 - Application version data (16 bytes)
short    nMajorVersion;
short    nMinorVersion;
short    nBugfixVersion;
short    nBuildVersion;
char     sUnused019[8];

// GROUP #19 - LTP protocol (14 bytes)
short    nLTPTType;
short    nLTPUsageOfDAC[ABF_WAVEFORMCOUNT];
short    nLTPPresynapticPulses[ABF_WAVEFORMCOUNT];
char     sUnused020[4];

// GROUP #20 - Digidata 132x Trigger out flag. (8 bytes)
short    nDD132xTriggerOut;
char     sUnused021[6];

// GROUP #21 - Epoch resistance (40 bytes)
char     sEpochResistanceSignalName[ABF_WAVEFORMCOUNT][ABF_ADCCOUNT];
short    nEpochResistanceState[ABF_WAVEFORMCOUNT];
char     sUnused022[16];

// GROUP #22 - Alternating episodic mode (58 bytes)
short    nAlternateDACOutputState;
short    nAlternateDigitalValue[ABF_EPOCHCOUNT];
short    nAlternateDigitalTrainValue[ABF_EPOCHCOUNT];
short    nAlternateDigitalOutputState;
char     sUnused023[14];

// GROUP #23 - Post-processing actions (210 bytes)
float    fPostProcessLowpassFilter[ABF_ADCCOUNT];
char     nPostProcessLowpassFilterType[ABF_ADCCOUNT];

// 6014 header bytes allocated + 130 header bytes not allocated
char     sUnused2048[130];

ABFFFileHeader();
}; // Size = 6144
// This structure is persisted, so the size MUST NOT CHANGE
STATIC_ASSERT(sizeof(ABFFFileHeader) == 6144);

inline ABFFFileHeader::ABFFFileHeader()
{
    // Set everything to 0.
    memset( this, 0, sizeof(ABFFFileHeader) );

    // Set critical parameters so we can determine the version.
    lFileSignature      = ABF_NATIVE_SIGNATURE;
    fFileVersionNumber  = ABF_CURRENT_VERSION;
    fHeaderVersionNumber = ABF_CURRENT_VERSION;
    lHeaderSize         = ABF_HEADER_SIZE;
}

//
// Scope descriptor format.
//
#define ABF_FACESIZE 32
struct ABFLogFont
{
    short nHeight;           // Height of the font in pixels.
    // short lWidth;           // use 0
    // short lEscapement;       // use 0
    // short lOrientation;      // use 0
    short nWeight;           // MSWindows font weight value.
    // char bItalic;           // use 0
    // char bUnderline;         // use 0
    // char bStrikeOut;         // use 0
    // char cCharSet;           // use ANSI_CHARSET (0)
    // char cOutPrecision;       // use OUT_TT_PRECIS
    // char cClipPrecision;      // use CLIP_DEFAULT_PRECIS
    // char cQuality;           // use PROOF_QUALITY
    char cPitchAndFamily;     // MSWindows pitch and family mask.
    char Unused[3];           // Unused space to maintain 4-byte packing.
    char szFaceName[ABF_FACESIZE]; // Face name of the font.
}; // Size = 40

struct ABFSignal
{

```

```

char    szName[ABF_ADCNAMELEN+2];        // ABF name length + '\0' + 1 for alignment.
short   nMxOffset;                       // Offset of the signal in the sampling sequence.
DWORD   rgbColor;                        // Pen color used to draw trace.
char     nPenWidth;                      // Pen width in pixels.
char     bDrawPoints;                   // TRUE = Draw disconnected points
char     bHidden;                       // TRUE = Hide the trace.
char     bFloatData;                    // TRUE = Floating point pseudo channel
float    fVertProportion;                // Relative proportion of client area to use
float    fDisplayGain;                   // Display gain of trace in UserUnits
float    fDisplayOffset;                 // Display offset of trace in UserUnits

// float    fUUTop;                      // Top of window in UserUnits
// float    fUUBottom;                   // Bottom of window in UserUnits
};    // Size = 34

////////////////////////////////////
//// WARNING WARNING WARNING WARNING WARNING WARNING WARNING WARNING WARNING ////
////////////////////////////////////
// The Following #defines appear to be largely unused in opur code base
// However there does exist a second set of #defines in AxScope32.h
// that REALLY defines what these bits in the header do.
// In particular it important to note that all 32 bits are in fact used internally
////////////////////////////////////
// Bit flags used in dwFlags field of ABFScopeConfig.
#define ABF_OVERLAPPED      0x00000001
#define ABF_DONTERASE       0x00000002
#define ABF_MONOCHROME      0x00000004
#define ABF_CLIPPING        0x00000008
#define ABF_HIDEHORZGRIDS   0x00000010
#define ABF_HIDEVERTGRIDS   0x00000020
#define ABF_FULLSCREEN      0x00000040
#define ABF_HIDEEXAXIS      0x00000080
#define ABF_HIDEYAXIS       0x00000100
#define ABF_HIDEXSCROLL     0x00000200
#define ABF_HIDEYSCROLL     0x00000400
#define ABF_HIDESIGNALNAME  0x00000800
#define ABF_ENABLEZOOM      0x00001000
#define ABF_XSPINFROMCENTER 0x00002000
#define ABF_HIDEEXSPINNER   0x00004000
#define ABF_LARGESPINNERS   0x00008000
#define ABF_PERSISTENCEMODE 0x00010000
#define ABF_CARDIACMODE     0x00020000
#define ABF_HIDETWIRLER     0x00040000
#define ABF_DISABLEUI       0x00080000
////////////////////////////////////
// #define ABF_INTERNALUSE  0xFFFF0000
// Do not add extra bit flags ^^^ here they are used internally
////////////////////////////////////
//// DANGER DANGER DANGER DANGER DANGER DANGER DANGER DANGER DANGER DANGER DANGER////
////////////////////////////////////

// Values for the wScopeMode field in ABFScopeConfig.
#define ABF_EPISODICMODE    0
#define ABF_CONTINUOUSMODE  1
// #define ABF_XYMODE        2

// Values for the nEraseStrategy field in ABFScopeConfig.
#define ABF_ERASE_EACHSWEEP  0
#define ABF_ERASE_EACHRUN    1
#define ABF_ERASE_EACHTRIAL  2
#define ABF_ERASE_DONTERASE  3

// Indexes into the rgbColor field of ABFScopeConfig.
#define ABF_BACKGROUND_COLOR 0
#define ABF_GRID_COLOR       1
#define ABF_THRESHOLD_COLOR  2
#define ABF_EVENTMARKER_COLOR 3
#define ABF_SEPARATOR_COLOR  4
#define ABF_AVERAGE_COLOR    5
#define ABF_OLD_DATA_COLOR    6
#define ABF_TEXT_COLOR        7
#define ABF_AXIS_COLOR        8
#define ABF_ACTIVE_AXIS_COLOR 9
#define ABF_LAST_COLOR        ABF_ACTIVE_AXIS_COLOR
#define ABF_SCOPE_COLORS      (ABF_LAST_COLOR+1)

// Extended colors for rgbColorEx field in ABFScopeConfig

```

```

#define ABF_STATISTICS_REGION0 0
#define ABF_STATISTICS_REGION1 1
#define ABF_STATISTICS_REGION2 2
#define ABF_STATISTICS_REGION3 3
#define ABF_STATISTICS_REGION4 4
#define ABF_STATISTICS_REGION5 5
#define ABF_STATISTICS_REGION6 6
#define ABF_STATISTICS_REGION7 7
#define ABF_BASELINE_REGION 8
#define ABF_STOREDSWEEPCOLOR 9
#define ABF_LASTCOLOR_EX ABF_STOREDSWEEPCOLOR
#define ABF_SCOPECOLORS_EX (ABF_LASTCOLOR+1)

// Values for the nDockState field in ABFScopeConfig
#define ABF_SCOPE_NOTDOCKED 0
#define ABF_SCOPE_DOCKED_TOP 1
#define ABF_SCOPE_DOCKED_LEFT 2
#define ABF_SCOPE_DOCKED_RIGHT 3
#define ABF_SCOPE_DOCKED_BOTTOM 4

struct ABFScopeConfig
{
    // Section 1 scope configurations
    DWORD dwFlags; // Flags that are meaningful to the scope.
    DWORD rgbColor[ABF_SCOPECOLORS]; // Colors for the components of the scope.
    float fDisplayStart; // Start of the display area in ms.
    float fDisplayEnd; // End of the display area in ms.
    WORD wScopeMode; // Mode that the scope is in.
    char bMaximized; // TRUE = Scope parent is maximized.
    char bMinimized; // TRUE = Scope parent is minimized.
    short xLeft; // Coordinate of the left edge.
    short yTop; // Coordinate of the top edge.
    short xRight; // Coordinate of the right edge.
    short yBottom; // Coordinate of the bottom edge.
    ABFLogFont LogFont; // Description of current font.
    ABFSignal TraceList[ABF_ADCCOUNT]; // List of traces in current use.
    short nYAxisWidth; // Width of the YAxis region.
    short nTraceCount; // Number of traces described in TraceList.
    short nEraseStrategy; // Erase strategy.
    short nDockState; // Docked position.
    // Size 656
    // * Do not insert any new members above this point! *
    // Section 2 scope configurations for file version 1.68.
    short nSizeofOldStructure; // Unused byte to determine the offset of the
version 2 data.
    DWORD rgbColorEx[ ABF_SCOPECOLORS_EX ]; // New color settings for stored sweep and cursors.
    short nAutoZeroState; // Status of the autozero selection.
    DWORD dwCursorsVisibleState; // Flag for visible status of cursors.
    DWORD dwCursorsLockedState; // Flag for enabled status of cursors.
    char sUnsigned[61];
    // Size 113
    ABFScopeConfig();
}; // Size = 769

inline ABFScopeConfig::ABFScopeConfig()
{
    // Set everything to 0.
    memset( this, 0, sizeof(ABFScopeConfig) );

    // Set critical parameters so we can determine the version.
    nSizeofOldStructure = 656;
}

//
// Definition of the ABF synch array structure
//

struct ABFSynch
{
    long lStart; // Start of the episode/event in fSynchTimeUnit units.
    long lLength; // Length of the episode/event in multiplexed samples.
}; // Size = 8

//
// Constants for nTagType in the ABFTag structure.
//
#define ABF_TIMETAG 0
#define ABF_COMMENTTAG 1
#define ABF_EXTERNALTAG 2

```

```

#define ABF_VOICETAG          3
#define ABF_NEWFILETAG        4

//
// Definition of the ABF Tag structure
//
struct ABFTag
{
    long    lTagTime;          // Time at which the tag was entered in fSynchTimeUnit units.
    char     sComment[ABF_TAGCOMMENTLEN]; // Optional tag comment.
    short    nTagType;         // Type of tag ABF_TIMETAG, ABF_COMMENTTAG, ABF_EXTERNALTAG or
ABF_VOICETAG.
    short    nVoiceTagNumber;  // If nTagType=ABF_VOICETAG, this is the number of this voice tag.
}; // Size = 64

// Comment inserted for externally acquired tags (expanded with spaces to ABF_TAGCOMMENTLEN).
#define ABF_EXTERNALTAGCOMMENT "<External>"
#define ABF_VOICETAGCOMMENT   "<Voice Tag>"

//
// Constants for nCompressionType in the ABFVoiceTagInfo structure.
//
#define ABF_COMPRESSION_NONE    0
#define ABF_COMPRESSION_PKWARE  1
// #define ABF_COMPRESSION_MPEG    2

//
// Definition of the ABFVoiceTagInfo structure.
//
struct ABFVoiceTagInfo
{
    long    lTagNumber;        // The tag number that corresponds to this VoiceTag
    long    lFileOffset;       // Offset to this tag within the VoiceTag block
    long    lUncompressedSize; // Size of the voice tag expanded.
    long    lCompressedSize;   // Compressed size of the tag.
    short    nCompressionType; // Compression method used.
    short    nSampleSize;      // Size of the samples acquired.
    long    lSamplesPerSecond; // Rate at which the sound was acquired.
    DWORD    dwCRC;            // CRC used to check data integrity.
    WORD     wChannels;        // Number of channels in the tag (usually 1).
    WORD     wUnused;          // Unused space.
}; // Size 32

//
// Constants for lParameterID in the ABFDelta structure.
//
// NOTE: If any changes are made to this list, the code in ABF_UpdateHeader must
//       be updated to include the new items.
#define ABF_DELTA_HOLDING0      0
#define ABF_DELTA_HOLDING1      1
#define ABF_DELTA_HOLDING2      2
#define ABF_DELTA_HOLDING3      3
#define ABF_DELTA_DIGITALOUTS    4
#define ABF_DELTA_THRESHOLD      5
#define ABF_DELTA_PRETRIGGER      6

// Because of lack of space, the Autosample Gain ID also contains the ADC number.
#define ABF_DELTA_AUTOSAMPLE_GAIN 100 // +ADC channel.

// Because of lack of space, the Signal Gain ID also contains the ADC number.
#define ABF_DELTA_SIGNAL_GAIN     200 // +ADC channel.

//
// Definition of the ABF Delta structure.
//
struct ABFDelta
{
    long    lDeltaTime;        // Time at which the parameter was changed in fSynchTimeUnit units.
    long    lParameterID;      // Identifier for the parameter changed
    union
    {
        long    lNewParamValue; // Depending on the value of lParameterID
        float    fNewParamValue; // this entry may be either a float or a long.
    };
}; // Size = 12

#ifdef RC_INVOKED
#pragma pack(pop) // return to default packing

```

```

#endif

//
// The size of the buffers to be passed to ABFH_GetWaveformVector
//
#define ABFH_MAXVECTORS      30

//
// Function prototypes for functions in ABFHADDR.C
//

void WINAPI ABFH_Initialize( ABFFileHeader *pFH );

void WINAPI ABFH_InitializeScopeConfig(const ABFFileHeader *pFH, ABFScopeConfig *pCfg);

BOOL WINAPI ABFH_CheckScopeConfig(ABFFileHeader *pFH, ABFScopeConfig *pCfg);

void WINAPI ABFH_GetADCDisplayRange( const ABFFileHeader *pFH, int nChannel,
                                     float *pfUUTop, float *pfUUBottom);

void WINAPI ABFH_GetADCToUUFactors( const ABFFileHeader *pFH, int nChannel,
                                     float *pfADCToUUFactor, float *pfADCToUUShift );
void WINAPI ABFH_ClipADCValue(const ABFFileHeader *pFH, int nChannel, float *pfUUVValue);

void WINAPI ABFH_GetDACtoUUFactors( const ABFFileHeader *pFH, int nChannel,
                                     float *pfDACtoUUFactor, float *pfDACtoUUShift );
void WINAPI ABFH_ClipDACValue(const ABFFileHeader *pFH, int nChannel, float *pfUUVValue);

BOOL WINAPI ABFH_GetMathValue(const ABFFileHeader *pFH, float fA, float fB, float *pfRval);
int  WINAPI ABFH_GetMathChannelName(char *pszName, UINT uNameLen);

BOOL WINAPI ABFH_ParamReader( HANDLE hFile, ABFFileHeader *pFH, int *pnError );
BOOL WINAPI ABFH_ParamReaderEx( HANDLE hFile, ABFFileHeader *pFH, int *pnError );
BOOL WINAPI ABFH_ParamWriter( HANDLE hFile, ABFFileHeader *pFH, int *pnError );

BOOL WINAPI ABFH_GetErrorText( int nError, char *pszBuffer, UINT nBufferSize );

// ABFHWAVE.CPP

// Constants for ABFH_GetEpochLimits
#define ABFH_FIRSTHOLDING  -1
#define ABFH_LASTHOLDING   ABF_EPOCHCOUNT

// Return the bounds of a given epoch in a given episode. Values returned are ZERO relative.
BOOL WINAPI ABFH_GetEpochLimits(const ABFFileHeader *pFH, int nADCChannel, DWORD dwEpisode,
                                int nEpoch, UINT *puEpochStart, UINT *puEpochEnd,
                                int *pnError);

BOOL WINAPI ABFH_GetEpochLimitsEx(const ABFFileHeader *pFH, int nADCChannel, UINT uDACChannel, DWORD
dwEpisode,
                                int nEpoch, UINT *puEpochStart, UINT *puEpochEnd,
                                int *pnError);

// Get the offset in the sampling sequence for the given physical channel.
BOOL WINAPI ABFH_GetChannelOffset( const ABFFileHeader *pFH, int nChannel, UINT *puChannelOffset );

// This function forms the de-multiplexed DAC output waveform for the
// particular channel in the pBuffer, in DAC UserUnits.
BOOL WINAPI ABFH_GetWaveform( const ABFFileHeader *pFH, int nADCChannel, DWORD dwEpisode,
                              float *pBuffer, int *pnError);

BOOL WINAPI ABFH_GetWaveformEx( const ABFFileHeader *pFH, UINT uDACChannel, DWORD dwEpisode,
                              float *pBuffer, int *pnError);

// This function forms the de-multiplexed Digital output waveform for the
// particular channel in the pdwBuffer, as a bit mask. Digital OUT 0 is in bit 0.
BOOL WINAPI ABFH_GetDigitalWaveform( const ABFFileHeader *pFH, int nChannel, DWORD dwEpisode,
                                     DWORD *pdwBuffer, int *pnError);

// Returns vector pairs for displaying a waveform made up of epochs.
BOOL WINAPI ABFH_GetWaveformVector(const ABFFileHeader *pFH, DWORD dwEpisode, UINT uStart,
                                   UINT uFinish, float *pfLevels, float *pfTimes,
                                   int *pnVectors, int *pnError);

// Returns vector pairs for displaying the digital outs.
BOOL WINAPI ABFH_GetDigitalWaveformVector(const ABFFileHeader *pFH, DWORD dwEpisode, UINT uStart,
                                           UINT uFinish, DWORD *pdwLevels, float *pfTimes,
                                           int *pnVectors, int *pnError);

// Calculates the timebase array for the file.

```



```

void WINAPI ABFH_GetTimebase(const ABFFileHeader *pFH, float fTimeOffset, float *pfBuffer, UINT
uBufferSize);
void WINAPI ABFH_GetTimebaseEx(const ABFFileHeader *pFH, double dTimeOffset, double *pdBuffer, UINT
uBufferSize);

// Get the duration of the first holding period.
UINT WINAPI ABFH_GetHoldingDuration(const ABFFileHeader *pFH);

// Checks whether the waveform varies from episode to episode.
BOOL WINAPI ABFH_IsConstantWaveform(const ABFFileHeader *pFH);

BOOL WINAPI ABFH_IsConstantWaveformEx(const ABFFileHeader *pFH, UINT uDACChannel);

// Checks that the sample intervals in the header are valid.
BOOL WINAPI ABFH_CheckSampleIntervals(const ABFFileHeader *pFH, float fClockResolution, int
*pnError);

// Gets the closest sample intervals higher and lower than the passed interval.
void WINAPI ABFH_GetClosestSampleIntervals(float fSampleInterval, float fClockResolution,
int nOperationMode, float fMinPeriod, float fMaxPeriod,
float *pfHigher, float *pfLower);

// Sets up the list for the spinner to drive the sampling interval through.
UINT WINAPI ABFH_SetupSamplingList(UINT uNumChannels, float fMinPeriod, float fMaxPeriod,
float *pfIntervalList, UINT uListEntries);

// Get the full sweep length given the length available to epochs or vice-versa.
int WINAPI ABFH_SweepLenFromUserLen(int nUserLength, int nNumChannels);
int WINAPI ABFH_UserLenFromSweepLen(int nSweepLength, int nNumChannels);

// Converts a display range to the equivalent gain and offset factors.
void WINAPI ABFH_GainOffsetToDisplayRange(const ABFFileHeader *pFH, int nChannel,
float fDisplayGain, float fDisplayOffset,
float *pfUUTop, float *pfUUBottom);

// Converts a display range to the equivalent gain and offset factors.
void WINAPI ABFH_DisplayRangeToGainOffset(const ABFFileHeader *pFH, int nChannel,
float fUUTop, float fUUBottom,
float *pfDisplayGain, float *pfDisplayOffset);

// Converts a time value to a synch time count or vice-versa.
void WINAPI ABFH_SynchCountToMS(const ABFFileHeader *pFH, UINT uCount, double *pdTimeMS);
UINT WINAPI ABFH_MSToSynchCount(const ABFFileHeader *pFH, double dTimeMS);

// Gets the point at which the sampling interval changes if split clock.
UINT WINAPI ABFH_GetClockChange(const ABFFileHeader *pFH);

// Gets the duration of the Waveform Episode (in us), allowing for split clock etc.
void WINAPI ABFH_GetEpisodeDuration(const ABFFileHeader *pFH, double *pdEpisodeDuration);

// Gets the duration of a P/N sequence (in us), including settling times.
void WINAPI ABFH_GetPNDuration(const ABFFileHeader *pFH, double *pdPNDuration);

void WINAPI ABFH_GetPNDurationEx(const ABFFileHeader *pFH, UINT uDAC, double *pdPNDuration);

// Gets the duration of a pre-sweep train in us.
void WINAPI ABFH_GetTrainDuration(const ABFFileHeader *pFH, double *pdTrainDuration);

void WINAPI ABFH_GetTrainDurationEx (const ABFFileHeader *pFH, UINT uDAC, double *pdTrainDuration);

// Gets the duration of a whole meta-episode (in us).
void WINAPI ABFH_GetMetaEpisodeDuration(const ABFFileHeader *pFH, double *pdMetaEpisodeDuration);

// Gets the start to start period for the episode in us.
void WINAPI ABFH_GetEpisodeStartToStart(const ABFFileHeader *pFH, double *pdEpisodeStartToStart);

// Checks that the user list contains valid entries for the protocol.
BOOL WINAPI ABFH_CheckUserList(const ABFFileHeader *pFH, int *pnError);

BOOL WINAPI ABFH_CheckUserListEx(const ABFFileHeader *pFH, UINT uListNum, int *pnError);

// Checks if the ABFFileHeader is a new (6k) or old (2k) header.
BOOL WINAPI ABFH_IsNewHeader(const ABFFileHeader *pFH);

// Demotes a 1.5 or 1.6 (or greater) ABF header to 1.5 version ABF header.
void WINAPI ABFH_DemoteHeader(ABFFileHeader *pOut, const ABFFileHeader *pIn );

// Promotes a 1.5 or 1.6 (or less) ABF header to a 1.6 ABF header.
void WINAPI ABFH_PromoteHeader(ABFFileHeader *pOut, const ABFFileHeader *pIn );

```

```

// Gets the first sample interval, expressed as a double.
double WINAPI ABFH_GetFirstSampleInterval( const ABFFileHeader *pFH );

// Gets the second sample interval expressed as a double.
double WINAPI ABFH_GetSecondSampleInterval( const ABFFileHeader *pFH );

// Counts the number of changing sweeps.
UINT WINAPI ABFH_GetNumberOfChangingSweeps( const ABFFileHeader *pFH );

// // Checks whether the digital output varies from episode to episode.
BOOL WINAPI ABFH_IsConstantDigitalOutput(const ABFFileHeader *pFH);

BOOL WINAPI ABFH_IsConstantDigitalOutputEx(const ABFFileHeader *pFH, UINT uDACChannel);

//
// Error return values that may be returned by the ABFH_xxx functions.
//

#define ABFH_FIRSTERRORNUMBER          2001
#define ABFH_EHEADERREAD                2001
#define ABFH_EHEADERWRITE               2002
#define ABFH_EINVALIDFILE               2003
#define ABFH_EUNKNOWNFILETYPE           2004
#define ABFH_CHANNELNOTSAMPLED          2005
#define ABFH_EPOCHNOTPRESENT             2006
#define ABFH_ENOWAVEFORM                 2007
#define ABFH_EDACFILEWAVEFORM           2008
#define ABFH_ENOMEMORY                   2009
#define ABFH_BADSAMPLEINTERVAL           2010
#define ABFH_BADSECONDSAMPLEINTERVAL     2011
#define ABFH_BADSAMPLEINTERVALS          2012
#define ABFH_ENOCONDITTRAINS             2013
#define ABFH_EMETADURATION               2014
#define ABFH_ECONDITNUMPULSES            2015
#define ABFH_ECONDITBASEDUR              2016
#define ABFH_ECONDITBASELEVEL            2017
#define ABFH_ECONDITPOSTTRAINDUR         2018
#define ABFH_ECONDITPOSTTRAINLEVEL       2019
#define ABFH_ESTART2START                 2020
#define ABFH_EINACTIVEHOLDING             2021
#define ABFH_EINVALIDCHARS               2022
#define ABFH_ENODIG                      2023
#define ABFH_EDIGHOLDLEVEL                2024
#define ABFH_ENOPNPULSES                  2025
#define ABFH_EPNNUMPULSES                 2026
#define ABFH_ENOEPOCH                     2027
#define ABFH_EEPOCHLEN                    2028
#define ABFH_EEPOCHINITLEVEL              2029
#define ABFH_EDIGLEVEL                    2030
#define ABFH_ECONDITSTEPDUR               2031
#define ABFH_ECONDITSTEPLEVEL             2032
#define ABFH_EINVALIDBINARYCHARS         2033
#define ABFH_EBADWAVEFORM                 2034

#ifdef __cplusplus
}
#endif

#endif /* INC_ABFHEADR_H */

```

## See Also:

[The ABF File I/O Functions](#)

[The ABF File I/O Functions by category](#)

[ABFFILES.H](#)

[ABFINFO.H](#)