



A signal analysis environment for  
MATLAB

Neuroscience Toolkit

Version 0.8 12<sup>th</sup> June 2008  
Author: Malcolm Lidieth  
King's College London

The sigTOOL Neuroscience Toolkit adds import functions for a number of commonly used commercial file formats to sigTOOL. It also provides analysis functions for spike train analysis and a link to Wave\_clus spike sorting software.

### Acknowledgements

Import of most of the formats below depends on the use of a redistributable Windows application extension available from the manufacturer's website. In addition, some use the MATLAB/Neuroshare API.

I am grateful to Professor Rodrigo Quian Quiroga of the University of Leicester, UK for permission to include an interface to the Wave-clus software.


File extension	Manufacturer	Software	Platforms
ABF <sup>1</sup>	<a href="#">Molecular Devices Inc</a> (Axon Instruments)	pClamp, AxoScope, ClampFit	Windows only
CFS <sup>2</sup>	<a href="#">Cambridge Electronic Design Ltd</a>	Signal	Windows only
MAP <sup>2</sup>	<a href="#">Alpha Omega</a>		Windows only
SMR/SON <sup>3</sup>	<a href="#">Cambridge Electronic Design Ltd</a>	Spike2 software	All
NEX <sup>3</sup>	<a href="#">Nex Technologies</a>	NeuroExplorer	Windows only
NEV <sup>3</sup>	<a href="#">Cyberkinetics Inc</a>		Windows Only
PLX <sup>3</sup>	<a href="#">Plexon Instruments</a>		Windows only

<sup>1</sup> Supported via the manufacturer's Windows application extension

<sup>2</sup> Supported via the manufacturer's Windows application extension and the Neuroshare/MATLAB API. See [www.neuroshare.org](http://www.neuroshare.org)

<sup>3</sup> Uses MATLAB native code.

The following analysis functions are added by the Neuroscience Toolkit to a Spike menu in sigTOOL data views:



These menus provide analysis functions for spike trains. They also allow a spike train to be converted to a rate channel that may be further analyzed using the standard Waveform menu analysis functions.

See the online help for further details.